



VWR INTERNATIONAL

ALL YOU NEED FOR SPECTROPHOTOMETRY

Life science applications

Applications in the environmental
and food industries

Universal photometers

Single beam devices

Double beam devices

Fluorometers

Turbidimeters for turbidity
measures

Cuvettes, accessories
and software

Quality assurance

MAKING SPECTROPHOTOMETRY *EASY!*

Introduction



Dear Customer,

The following questions and answers are designed to help you decide which colorimeter, filter or spectrophotometer is best suited to your requirements. Please remember that a matrix like this cannot always replace the experience and competence that an advisor can offer. You can therefore contact our field specialists if you have any further questions.

Photometry is an important part of modern instrumental analysis. It is now an established measurement and testing method in pharmaceutical quality control, in water and environmental laboratories and in quality assurance in the food industry. Photometers have as many individual features as they do areas of application. We offer an extensive range of devices so we can meet all your requirements in educational institutions, industry and in public and private research laboratories.

Our portfolio starts with simple colorimeters and goes right through to spectrophotometers and the related complete test kits for rapid and precise analysis. A summary table gives you a rapid overview of our product range. We have compiled a list of questions to guide your decision so you can choose the most appropriate instrument for your specific application. Our range is completed by standard solutions and test systems for your quality management, tips for accessories and consumables, which allows us to help you ensure high quality standards. Our field product specialists will also be happy to demonstrate the instruments on site for you.

Introduction

| | |
|---|-----|
| Choosing the most appropriate photometer | 3 |
| Instrument overview | 4-7 |

Photometers for life science applications

| | |
|--|----|
| Spectrophotometer GeneQuant 1300/100..... | 8 |
| Cell density meter Ultrospec 10 | 8 |
| Photometers Ultrospec 2100 pro/3100 pro..... | 9 |
| Photometers Novaspec Plus..... | 9 |
| Photometer GENESYS™ 10S BIO UV/VIS | 10 |
| Biophotometer plus..... | 11 |
| Photometers Novaspec III | 11 |
| Spectrophotometer Nanovue™ Plus..... | 12 |
| Microplate reader Multiskan® FC..... | 13 |
| Microplate reader Multiskan® GO | 14 |

Photometers for environmental and food industry applications

| | |
|--|----|
| Colorimeter Spectroquant® Picco | 15 |
| Colorimeter Spectroquant® Multy | 16 |
| Photometers Spectroquant® | |
| NOVA 30A/NOVA 60/60A..... | 17 |
| Spectrophotometers | |
| Spectroquant® Pharo 100/Pharo 300..... | 18 |
| Spectrophotometer AquaMate™ plus..... | 19 |

Universal photometers - single beam devices

| | |
|---|----|
| Spectrophotometers VWR Collection V-1200 ... | 20 |
| Spectrophotometer GENESYS™ 20 VIS | 21 |
| Spectrophotometers | |
| VWR Collection UV-1600 PC..... | 22 |
| Spectrophotometers | |
| VWR Collection UV-3100 PC..... | 23 |
| Spectrophotometers | |
| GENESYS™ 10S VIS UV/VIS..... | 24 |
| GENESYS™ photometers software | 25 |
| Photometers Ultrospec 5300 pro/6300 pro..... | 25 |
| Spectrophotometers HITACHI U-5100..... | 26 |
| Spectrophotometers | |
| UviLine 9100/9400 UV/VIS | 27 |
| Photometers single beam devices | |
| JENWAY 7300 series | 28 |
| Spectrophotometer split beam 6700 series..... | 29 |

Universal photometers - double beam devices

| | |
|--------------------------------|----|
| Spectrophotometers double beam | |
| UV/VIS U-2900/U-2910 | 30 |
| Spectrophotometers double beam | |
| UV/VIS U-3900/U-300H..... | 31 |

Fluorometers

| | |
|--|----|
| Fluorescence spectrophotometer F-2700..... | 32 |
| Fluorescence spectrophotometer F-7000..... | 33 |

Turbidimeters for measuring turbidity

| | |
|--|-------|
| Turbidimeters Turbiquant® 1100 IR/T, 1500 IR/T, 3000 IR/T | 34-35 |
|--|-------|

Cuvettes

| | |
|--|-------|
| Absorption measurement macro standard cuvettes..... | 36 |
| Semi-micro, micro cuvettes | 37 |
| Disposable cuvettes | 38-39 |
| Accessories, cell stands, cleaning agents | 39 |

Quality assurance

| | |
|--|----|
| Secondary standards for calibrating spectrophotometers | 40 |
| AQS-1 mode for monitoring photometers, UV/VIS standards Certipur® | 41 |
| AQS-2 mode standard solutions | |
| Spectroquant® CombiCheck..... | 42 |
| AQS-3 mode standard solutions | |
| Spectroquant® | 43 |

| | |
|---|---|
| What field of application do you operate in? | BioTech/life science: GENESYS™ 10S BIO, Biophotometer plus, Novaspec series, Ultrospec series, GeneQuant |
| | Water/wastewater analysis: Spectroquant® NOVA, Pharo, Picco and Multy are suitable for Spectroquant® test kits from Merck, whilst the AquaMate™ plus is suitable as an independent device |
| | Pharmaceutical (regulated): Model 6715, U-2900/2910, U-3900/3900H, Ultrospec 5300/6300 |
| Do you only want to take measurements in the visual range (325 - 1100 nm) or in the UV range (190 - 325 nm) as well? | General material research (colour measurement): U-3900/3900H |
| | Visual measuring range: VWR V-1200, GENESYS™ 20, GENESYS™ 10S VIS, AquaMate™ plus VIS, Model 6700, Models 7300/7305, UviLine 9100 VIS |
| Which sample volumes do you have available? | UV and visual measuring range: All other UV/VIS spectrophotometers |
| | Macro/micro (>300 µl filling volumes): All devices are suitable |
| Do you want to save more than one method or measurement data set in the device or on an external medium (not a PC)? | Ultra-micro (20 - 200 µl filling volumes): Special plastic single use cuvettes with approx. 50 µl filling volume must be handled separately - the cuvette must be securely located in the cell holder, a pin hole is used ideally to deliver reproducible results; do not use an automated cell changer as the reproducibility of the measuring point is not given |
| | USB memory sticks: UviLine 9100, 9400, U-2900, AquaMate™ plus, all GENESYS™ 10S models, Models 7310/7315 |
| Do you want to print the data/results? | SD card: Models 6700/6705/6715 |
| | Internal printer: GENESYS™ 20/10S, 10S BIO, 6700/6705/6715, 7300/7305/7310/7315, GeneQuant 100/1300, NanoVue can be ordered with an integrated thermal transfer printer |
| Do you want to measure the sample directly from/in a sample vessel (and not transfer the samples to a cuvette)? | External printer port: All devices taking varying printer specifications into consideration |
| | Sipper (a pump is used to suck the sample into a pre-installed flow through cell): A range of manual or automatic sipper systems are available for many of the units, please check the accessory tables |
| Do you want to measure more than one sample at the same time (possibly with zero compensation)? | Glass fibre measuring head connection (the measuring sensor is immersed directly into the sample solution outside the photometer): GENESYS™ 10S/10S BIO |
| | Automatic multi-cell changer: VWR UV-1600PC, UV-3100PC, Models 67 and 73 series, UviLine 9100/9400, GENESYS™ 10S/10S BIO, AquaMate™ plus, U-5100 |
| | Automatic multi-cell changer with Peltier cooling: U-2900/2910, U-3900 |
| Do you want to work with standards (create a calibration curve)? | Autosampler: AquaMate™ plus, U-2900/2910, U-3900 |
| Are your samples highly diluted? | All devices |
| Do you want to constantly measure against a reference sample? | All devices: Special cell holder 20 - 50/100 mm cuvettes |
| Would you like to operate a stand-alone instrument by a PC? | Double beam devices: U-2900/2910, U-3900 |
| Do you want to carry out fluorescence measurements? | VWR UV-1600PC, UV-3100PC, GENESYS™ 20/10S/10S BIO, AquaMate™ plus, U-5100, U-2900, 73 series |
| | Fluorescence spectrophotometers F-2700, F-7000 |

Introduction

● Photometers for life science applications

| Description | Light source (lamps) | Optical system | Wavelength range (nm) (2) | Spectral bandwidth (nm) | Photometric range (A) |
|---|----------------------|-------------------------------------|--|-------------------------|-----------------------|
| GeneQuant 1300/100 | Xenon | Single beam with internal reference | 190 - 1100 | 5 | -0,3 to +2,5 |
| Ultrospec 10 | LED | Single beam | 600 | 40 | -0,3 to 1,99 |
| Ultrospec 2100 pro, Ultrospec 3100 pro | Xenon | Single beam with internal reference | 190 - 1100 | <3 | -3,0 to +3,0 |
| Novaspec plus, Novaspec III | Tungsten | Single beam | 330 - 800 | 7 | -0,3 to +2,5 |
| Ultrospec 5300 pro, Ultrospec 6300 pro | Tungsten, deuterium | Single beam with internal reference | 190 - 1100 | 1 | -3,0 to +3,0 |
| GENESYS™ 10S UV, GENESYS™ 10S BIO | Xenon | Split beam with internal reference | 190 - 1100 | 1,8 | -0,1 to +3,0 |
| BioPhotometer plus | Xenon | Single beam with internal reference | 230, 260, 280, 340, 405, 490, 550, 595, 650 | 5 - 7 | 0 to 3 |
| NanoVue™ Plus | Xenon | Single beam with internal reference | 190 - 1100 | 5 | N/A |
| Multiskan® FC | Tungsten | Single beam | 340, 405, 414, 492, 540, 620 | 3 - 9 (½ bandwidth) | 0 to 4 |
| Multiskan® GO | Xenon | Single beam | 200 - 1000 | <2,5 | 0 to 2,5 |

● Photometers for applications in the environmental and food industries

| Description | Light source (lamps) | Optical system | Wavelength range (nm) (2) | Spectral bandwidth (nm) | Photometric range (A) |
|---|----------------------|-------------------------------------|---|--|-----------------------|
| Spectroquant® Picco Cl ₂ / O ₃ / ClO ₂ / CyA / pH | LED | - | (1) | Single parameter devices for the COD, fluoride, ammonium, pH, | |
| Colorimeter Spectroquant® Multy | 6x LED | Single beam with internal reference | 430, 530, 560, 580, 610, 660 | | |
| Spectroquant® NOVA 30A | Tungsten | Single beam with internal reference | 340, 445, 525, 550, 605, 690 | | -0,3 to +3,2 |
| Spectroquant® NOVA 60A | Tungsten | Single beam with internal reference | 340, 410, 445, 500, 525, 550, 565, 605, 620, 665, 690, 820 | | -0,3 to +3,2 |
| Spectroquant® Pharo 100 | Tungsten | Single beam, stabilised | 320 - 1100 | 4 | ±3,3 |
| Spectroquant® Pharo 300 | Xenon | Single beam, stabilised | 190 - 1100 | 4 | ±3,3 |
| AquaMate™ plus VIS | Tungsten | Single beam | 315 - 1100 | 2 | -0,3 to +3,0 |
| AquaMate™ plus UV/VIS | Tungsten, deuterium | Single beam | 190 - 1100 | 2 | -0,3 to +3,0 |

(*) PC required (no stand-alone functionality)

(1) Details for wavelengths and product codes see product page

(2) If single values mentioned main units are filter photometers

Device overview

| Interfaces | Memory internal/external | Specialities | Cat. No. |
|---|---------------------------------|--|---|
| USB | Yes/- | Predefined life science methods nucleic acid, protein assay, cell culture Cy™ Dye (only Genequant 1300) | 28-9182-04, 28-9182-13 |
| N/A | Yes/- | OD600 measurements (cell density) | 80-2116-30 |
| RS232, parallel | -/- | Predefined life science methods, cDNA (only Ultrospec-3100 pro) | 80-2112-21, 80-2112-31 |
| RS232 | -/- | Novaspec III basic unit for education, Novaspec plus basic unit for enzyme methods | 80-2117-50, 80-2118-00 |
| RS232, centronics (Ultrospec 6300 only) | -/- | FDA compliant software SWIFT II available (optional) | 80-2117-70, 80-2117-60 |
| USB | | Comprehensive software package, integral printer (optional), predefined life science methods (G10S BIO only) | 634-0595, 634-0596, 634-0591, 634-0592 |
| RS-232C | 100 results, 32 methods | 32 methods, 9 freely programmable, very small footprint | 634-4033, EU plug 634-0590, UK plug |
| Bluetooth®, SD card (optional) | 81 methods, SD card | Micro volume sample (min. 0,5 µl at 0,2 mm path length) | 28-9569-65 |
| USB | 99 protocols, results/USB stick | Ease of use, stand-alone mode | 736-0355, 736-0356 |
| USB | Yes/USB stick | Stand-alone mode, very fast, cuvette measurements | 736-0559, 736-0560 |

| Interfaces | Memory internal/external | Specialities | Cat. No. |
|---|----------------------------|---|-------------------------------|
| determination of chlorine, ozone, chlorine dioxide, cyanuric acid, pH – nitrate, nitrogen, o-phosphate/total phosphorous - see product page | | | (1) |
| RS232 for printer or PC | 1000 data sets | Multi-parameter device for 120 Spectroquant® cell tests | 1.73630.0001 |
| RS232C serial | 500 measured values/PC | Barcode reader for Spectroquant® cell tests | 1.09748.0001 |
| RS232C serial | 1000 measured values/PC | Barcode reader for Spectroquant® cell tests, 50 free programmable methods, spectral data and kinetics; NOVA 60A includes rechargeable batteries | 1.09751.0001, 1.09752.0001 |
| RS232, USB-A, USB-B | 1000 results/USB stick, PC | Barcode reader for Spectroquant® cell tests, freely programmable methods, spectral data and kinetics | 1.00706.001 |
| RS232, USB-A, USB-B | 1000 results/USB stick, PC | Barcode reader for Spectroquant® cell tests, freely programmable methods, spectral data and kinetics | 1.00707.0001 |
| RS232C, USB | Yes/USB stick | Comprehensive software package (optional) Uvcalc™ for flexible calculations | 705-0951 |
| RS232C, USB | Yes/USB stick | Comprehensive software package (optional) Uvcalc™ for flexible calculations | 705-0950 |

All units with PC remote control software can save methods and results on PC within remote software.

Introduction

● Universal photometers - single beam devices

| Description | Light source (lamps) | Optical system | Wavelength range (nm) (2) | Spectral bandwidth (nm) | Photometric range (A) |
|------------------|----------------------|---|---------------------------|-------------------------|-----------------------|
| VWR V-1200 | Tungsten | Single beam | 325 - 1000 | 4 | -0,3 to +3,0 |
| GENESYS™ 20 | Tungsten | Single beam | 325 - 1100 | 8 | -0,1 to +2,5 |
| VWR UV-1600PC | Tungsten, deuterium | Single beam | 190 - 1100 | 4 | -0,3 to +3,0 |
| VWR UV-3100PC | Tungsten, deuterium | Single beam | 190 - 1100 | 2 | -0,3 to +3,0 |
| GENESYS™ 10S VIS | Tungsten | Single beam | 325 - 1100 | 5 | -0,1 to +3,0 |
| U-5100 | Tungsten, deuterium | Ratio beam with internal reference | 190 - 1100 | 5 | -0,3 to +3,0 |
| UviLine 9100 | Tungsten | Single beam | 320 - 1000 | 4 | -3,3 to +3,3 |
| UviLine 9400 | Xenon | Single beam | 198 - 1000 | 4 | -3,3 to +3,3 |
| Model 7300 | Tungsten | Single beam | 320 - 1000 | 5 | -0,3 to +2,5 |
| Model 7305 | Xenon | Single beam | 198 - 1000 | 5 | -0,3 to +2,5 |
| Model 7310 | Tungsten | Single beam | 320 - 1000 | 5 | -0,3 to +2,5 |
| Model 7315 | Xenon | Single beam | 1908- 1000 | 5 | -0,3 to +2,5 |
| Model 6700 | Tungsten | Single split beam with internal reference | 320 - 1100 | 4 | -0,3 to +3,0 |
| Model 6705 | Xenon | Single split beam with internal reference | 190 - 1100 | 4 | -0,3 to +3,0 |
| Model 6715 | Xenon | Single split beam with internal reference | 190 - 1100 | 1,5 | -0,3 to +3,0 |

● Universal photometers - double beam devices

| Description | Light source (lamps) | Optical system | Wavelength range (nm) (2) | Spectral bandwidth (nm) | Photometric range (A) |
|--------------------|----------------------|----------------|---------------------------|-------------------------|-----------------------|
| U-2900, U-2910 (*) | Tungsten, deuterium | Double beam | 190 - 1100 | 1,5 | -0,3 to +3,0 |
| U-3900 (*) | Tungsten, deuterium | Double beam | 190 - 900 | 0,1 - 5 | -5,5 to +5,5 |
| U-3900H (*) | Tungsten, deuterium | Double beam | 190 - 900 | 0,1 - 5 | -3,8 to +3,8 |

● Fluorimeters

| Description | Light source (lamps) | Optical system | Wavelength range (nm) (2) | Spectral bandwidth (nm) | Photometric range (A) |
|-------------|----------------------|----------------|---------------------------|-------------------------|-----------------------|
| F-2700 | Xenon | | | | |
| F-7000 | Xenon | | | | |

(*) PC required (no stand-alone functionality)

(1) Details for wavelengths and product codes see product page

(2) If single values mentioned main units are filter photometers

Device overview

| Interfaces | Memory internal/external | Specialities | Cat. No. |
|----------------------------------|-------------------------------------|--|--------------------|
| RS232, USB | 200 results, 200 curves | | 634-6000 |
| RS232, parallel | | Integral printer (optional) | 634-1021, 634-1022 |
| RS232, USB | 200 results, 200 curves | Software for PC control included | 634-6001 |
| RS232, USB | 200 results, 200 curves | Full scan capabilities; software for PC control included | 634-6002 |
| USB | | Comprehensive software package, integral printer (optional) | 634-0593, 634-0594 |
| USB, parallel | 20 methods, 50 results | Comprehensive software package (optional) | 634-0757 |
| RS232, USB | 100 methods, 1000 results/USB stick | Optional extra keyboard, functional keys, admin./user system | 634-5000 |
| RS232, USB | 100 methods, 1000 results/USB stick | Optional extra keyboard, functional keys, admin./user system | 634-5001 |
| Analogue, RS232 | | Integral printer (optional), admin./user system | 664-0067 |
| USB, RS232, analogue | USB stick | Integral printer (optional), admin./user system | 664-0066 |
| Analogue, RS232 | | Integral printer (optional), admin./user system | 664-0065 |
| USB, RS232, analogue | USB stick | Integral printer (optional), admin./user system | 664-0064 |
| USB, parallel, analogue USB, USB | >1000 methods and results/SD card | Colour touch screen, Qwheel™, admin./user system | 634-0583 |
| USB, parallel, analogue USB, USB | >1000 methods and results/SD card | Colour touch screen, Qwheel™, admin./user system | 634-0584 |
| USB, parallel, analogue USB, USB | >1000 methods and results/SD card | Colour touch screen, Qwheel™, admin./user system | 634-0585 |

| Interfaces | Memory internal/external | Specialities | Cat. No. |
|--|------------------------------------|--|--------------------|
| RS232 to PC (U-2900/U-2910) parallel, USB (U-2900) | Yes/USB memory stick (U-2900 only) | Comprehensive software package in line with Ph. Eur. (optional) | 634-0718, 634-0719 |
| USB | | Comprehensive software package in line with Ph. Eur. (optional) | 634-0756 |
| USB | | Comprehensive software package in line with Ph. Eur. (optional), U-3900H: Double monochromator | 634-0754 |

| Interfaces | Memory internal/external | Specialities | Cat. No. |
|------------|--------------------------|--------------|----------|
| | | | 634-0067 |
| | | | 634-0743 |

All units with PC remote control software can save methods and results on PC within remote software.

Photometers for life science applications

Spectrophotometer GeneQuant 1300 / 100



Ordering information

| Description | Cat. No. |
|------------------------------------|------------|
| GeneQuant 1300 Classic | 28-9182-13 |
| GeneQuant 1300 with printer | 28-9182-14 |
| GeneQuant 1300 with Bluetooth® | 28-9182-15 |
| GeneQuant 100 Classic | 28-9182-04 |
| GeneQuant 100 Classic with printer | 28-9182-05 |

Compact, convenient and flexible.

The GeneQuant 1300 spectrophotometer is a compact, convenient and flexible instrument pre-programmed with a range of methods for the quantitation of proteins, nucleic acids, and bacterial cell cultures.

Absorbance and concentration measurements can be made at any wavelength, providing flexibility for future applications.

GeneQuant 100 is also available, with built-in applications for nucleic acids, proteins and cell density measurements.

- Built-in applications for cell density, nucleic acids, proteins, and enzyme kinetics
- cDNA application software for measuring incorporation of Cy™ 3 and Cy™ 5, or other dyes into purified microarray hybridisation probes and PCR products
- Wavelength range of 190 to 1100 nm, with wavelength scan from 200 to 900 nm performed in less than 5 s with zoom facility, peak identification, and on peak confirmation
- Visual inspection of a nucleic acid scan can identify the presence of impurities, especially useful with RNA samples
- Optional integrated printer
- USB cable and print via computer software supplied as standard - wireless Bluetooth® accessory (optional)
- Can be used with a wide variety of cuvettes or capillaries for sample volumes from 3 µl to 2 ml
- Three year lamp warranty

Cell density meter Ultrospec 10

A small, portable and dedicated cell density OD600 measurement instrument.

The Ultrospec 10 measures the density of cells (*E. coli*, other bacteria, and yeast*) in suspension at 600 nm. This battery powered, hand held device can be easily used in incubation cabinets, under anaerobic conditions and in other areas of the laboratory where cells are cultured. Rechargeable batteries provide nearly one month of cord-free use.

* not recommended for use with mammalian cell lines

- Battery operated for use in cell culture facilities
- Easy-to-use, clean and sterilise
- Download results to a PC or printer for convenient data storage and retrieval
- Available in four colours



Ordering information

| Description | Cat. No. |
|--|------------|
| Ultrospec 10 cell density meter, Classic | 80-2116-30 |
| Yellow | 80-2116-31 |
| Plum | 80-2116-32 |
| Apple | 80-2116-33 |

Photometers for life science applications

Photometer Ultrospec 2100 / 3100 pro



The Ultrospec 2100 pro is a simple to use instrument for nucleic acid measurements and enzyme kinetics.

It provides stored routines for nucleic acid quantitation and a standard curve routine for protein determination, in addition to measuring absorbance and concentration.

An 8 position sample changer is supplied as standard.

Accessories such as a temperature control unit option are also available.

The instrument can be upgraded for more sophisticated applications, as well as data storage, with SWIFT II software and a PC.

Ordering information

| Description | Cat. No. |
|-----------------------------|------------|
| Ultrospec 2100 pro, Classic | 80-2112-21 |
| Yellow | 80-2112-22 |
| Plum | 80-2112-27 |
| Apple | 80-2112-28 |
| Ultrospec 3100 pro, Classic | 80-2112-31 |
| Yellow | 80-2112-32 |
| Plum | 80-2112-33 |
| Apple | 80-2112-34 |

- Pre-set routines for DNA, RNA, oligonucleotide quantitation, and purity checks
- Combines basic measurement with graphic modes: Wavelength scans, absorbance changes with time, reaction rate determinations, and standard curves can be displayed and printed
- Stores up to 18 user-defined methods
- Includes serial and parallel outputs as standard (outputs to a range of printers)
- GLP self-test diagnostics
- IQ/OQ certification test plans
- Available in a choice of colours
- The Ultrospec 3100 is also available with all the benefits of the Ultrospec 2100, plus high resolution display, ability to store up to 50 user-defined methods and a cDNA application for measuring Cy[™] 3, Cy[™] 5 and fluorescein incorporation into probes and PCR products.

Photometers Novaspec Plus



A good value instrument for general laboratory use.

It includes stored methods for protein quantitation and enzyme kinetics, plus the basic modes of absorbance, transmittance, OD600, and concentration.

- Kinetics for enzyme studies
- Stored protein methods for Bradford, BCA, Biuret, and Lowry
- Bacterial cell culture measurement at OD600
- "Flash Scan" diode array, for rapid wavelength scans
- Graphic display of wavelength scans, kinetic assays (including slope calculation for rate/activity studies) and standard curves
- Option for temperature control (factory fitted)
- Up to 99 stored methods

Ordering information

| Description | Cat. No. |
|--|------------|
| Novaspec Plus | 80-2117-50 |
| Novaspec Plus with heated cell holder includes: Grafico PC utility software, serial cable and dust cover | 80-2117-51 |

Photometers for life science applications

UV/VIS photometer GENESYS™ 10S BIO



Delivery information:

Supplied with a 6 position cell changer, single cell holder, spare fuses, USB memory device, USB cable, dust cover and operator manual.

For accessories and software see GENESYS™ 10S VIS UV/VIS on page 24

GENESYS™ 10S BIO from Thermo Scientific

These upgraded models are compact, robust, economical, easy-to-use, scanning, split beam UV/VIS instruments with improved performance.

This instrument is versatile and includes intuitive and powerful internal software with pre-configured methods that can be edited and stored to a readily accessible customised assay.

A wide variety of accessories are available for the wide ranging requirements of life science applications (please ask for details).

- USB connectivity
- 1,8 nm bandwidth for improved nucleic acid measurement (ratio and concentration)
- Direct protein measurement at 280 and 205 nm
- Standard measurement of protein concentrations using Coomassie/Bradford, Lowry, modified Lowry, BCA, Pierce Micro-BCA and 660 Protein assays
- 6 cell changer for automated multi-sample experiments

Technical specifications

| | |
|--------------------------|--|
| Lamp source, lifetime | Xenon flash, 5 years typically |
| Optical system | Split beam with reference detector |
| Range (nm) | 190 - 1100 |
| Spectral bandwidth (nm) | 1,8 |
| Wavelength accuracy (nm) | ±1,0 |
| Scan speed (nm/min) | 10 - 4200 |
| Photometric range (A) | Up to 3,5 A at 260 nm |
| Photometric accuracy | ±0,005 A at 1,0 A 0,010 A K ₂ Cr ₂ O ₇ |
| Stray light | <0,08% T at 220, 340 nm (NaI, NaNO ₃) <1,0% 198 nm (KCl) |
| Display | Graphic 320 x 240 pixel backlit LCD (91 x 71 mm) |
| Standard cell holder | Integral 6 position cell changer, single cell holder |
| Keyboard | Membrane keypad |
| Printer | Optional, 40 column, internal, graphic |
| Interfaces | USB type A port for USB memory device (front panel) USB type B port for optional PC (rear panel) USB type A port for external printer (rear panel) |
| Power (V) | 100 - 240, selected automatically |
| W x D x H (mm) | 300 x 400 x 250 |
| Weight (kg) | 8,6 |

Ordering information

| Description | Cat. No. |
|---|----------|
| GENESYS™ 10S UV/VIS BIO, power leads with EU and UK plug | 634-0591 |
| GENESYS™ 10S UV/VIS BIO with internal printer and power leads with EU and UK plug | 634-0592 |

Photometers for life science applications

BioPhotometer plus



Compact and manageable device from Eppendorf®

This photometer is suitable for the rapid and reliable analysis of nucleic acids, proteins, cell density, dye, assays and general absorption measurements.

- Easy-to-use with a measurement time of approx. 2 seconds
- Small but compact device thanks to the robust metal case
- Choice of background compensation at 320 nm for UV measurements
- Direct calculation of all results taking dilutions into account
- 9 wavelengths for 32 methods, alterable

Ordering information

| Description | Cat. No. |
|-----------------------------|----------|
| BioPhotometer plus, EU plug | 634-4033 |
| BioPhotometer plus, UK plug | 634-0590 |

Accessories

| Description | Pk | Cat. No. |
|---|----|----------|
| DPU 414 thermal printer, including power supply and printer cable | 1 | 713-1001 |
| Thermal paper for printer | 5 | 700-5118 |
| Secondary UV/VIS test filter | 1 | 732-6015 |

Technical specifications

| | |
|-----------------------|---|
| Light source | Xenon flash |
| Optical system | Single beam optics with reference beam |
| Wavelength range (nm) | 230, 260, 280, 340, 405, 490, 550, 595, 650 |
| Spectral bandwidth | 5 nm at 230 - 340 nm, 7 nm at 405 - 650 nm |
| Wavelength accuracy | ±1 nm at 230 - 280 nm, ±2 nm at 550 - 650 nm |
| Photometric range | 0 - 3 A; 0 - 2 A at 340 nm, 0 - 2 A at Dye 550/650 nm |
| Display | Graphic LCD, illuminated (33x60 mm) |
| Keyboard | 19 membrane keys |
| Test storage | 32 pre-programmed (alterable) |
| Measurement methods | DNA, RNA, Protein, OD600, Dye, Assay, Absorption |
| Interfaces | RS232C serial for PC and printer |
| Power (V) | 100 - 240 |
| W x D x H (mm) | 200 x 320 x 100 |
| Weight (kg) | 3 |
| Warranty | 2 years |
| Accessories included | Dust cover, user manual, 8 cell UVette |

Photometers Novaspec III

A lightweight instrument for measuring absorbance, % transmission, concentration and rate.

It is so easy-to-use. Just select wavelength, set reference and measure sample, making it the instrument of choice for teaching laboratories.

- Large display with easy-to-read characters
- Absorbance, % transmission, concentration, factor and rate
- Ability to measure two wavelengths simultaneously for kinetics assays
- Self-test on start-up
- Analogue output for connection to chart recorder
- 3 student-friendly experiments plus a tutorial on UV/VIS spectrophotometry are included



Ordering information

| Description | Cat. No. |
|-------------------------------------|------------|
| Novaspec III | 80-2118-00 |
| Test tube adapters (10, 12 & 16 mm) | 80-2117-47 |

Photometers for life science applications

Spectrophotometer NanoVue™ Plus



Measurement of low volumes made easy

The NanoVue™ Plus is an easy-to-use and reliable instrument for the measurement of nucleic acid and protein samples. Samples of 0,5 to 2 µl can be pipetted directly onto a novel gold sample plate for measurement and then simply recovered using a pipette. If sample recovery is not required, the sample plate with hydrophobic coating can be quickly and easily wiped clean.

- Ideally suited to the life scientist where the sample is limited and speed and convenience of analysis is key
- Eliminates the need for cuvettes, capillaries or other sample devices - just drop and read
- Pre-defined methods for nucleic acid quantitation, including concentration, purity and theoretical T_m
- Pre-defined methods for protein quantitation, including direct UV, Bradford, BCA, Biuret and Lowry, with the ability to run up to 27 standards



Ordering information

| Description | Cat. No. |
|-------------------------------|------------|
| NanoVue™ Plus with SD card | 28-9569-60 |
| NanoVue™ Plus | 28-9569-65 |
| NanoVue™ Plus with printer | 28-9569-66 |
| NanoVue™ Plus with Bluetooth® | 28-9569-67 |

Supplied complete with software for Print Via Computer (PVC) and software Cy™ Dye and calibration fluid.

Technical specifications

| | |
|-------------------------------|--|
| Lamp source | Long life stabilised Xenon |
| Detector | Twin CCD array |
| Range (nm) | 200 - 1100 (scanning 200 - 950) |
| Spectral bandwidth (nm) | 5 |
| Wavelength accuracy (nm) | ±2 |
| Wavelength repeatability (nm) | ±0,5 |
| Photometric range (A) | 0 to 125 (10 mm path length equivalence) |
| Photometric accuracy | Maximum ±1% at 259 nm at 0,7 to 0,8 A using uracil |
| Method memory | 90 |
| Display | Backlit graphic LCD |
| Interfaces | USB port, optional Bluetooth® |
| Power | 100 to 240 VAC ±10%, 50/60 Hz, 50 VA |
| W x D x H (mm) | 260 x 390 x 100 |
| Weight (kg) | <4,5 |

Photometers for life science applications

Microplate reader Multiskan® FC Thermo Scientific



The Multiskan® FC, based on 30 years experience in the area of microplate measurement.

It is a compact, reliable and robust instrument capable of reading 96 and 384 well plates for a wide variety of applications.

It provides fast and accurate measurement with excellent linearity. One model in the range is fitted with an incubator providing temperature control up to 50 °C and the capacity to read 384 well plates. It takes only 6 seconds to measure a 96 well plate and 12 seconds to measure a 384 well plate. The spectral range of 340 to 850 nm makes the Multiskan® FC suited to applications from enzyme kinetics to Lowry assays.

- USB stick transfer of stored data with up to 100 assays in internal memory
- Use either as a stand-alone instrument using internal software or under PC control with SkanIt® software
- Small footprint saves laboratory space and facilitates automated procedures
- Reliable day-to-day and year-on-year performance
- Robotic compatibility for high throughput environments

Ordering information

| Description | Cat. No. |
|--|----------|
| Multiskan® FC with integral shaker | 736-0355 |
| Multiskan® FC with integral shaker and incubator | 736-0356 |

Supplied complete with 3 standard filters - (405, 450, 620 nm) and SkanIt® software

Technical specifications

| Model | Multiskan® FC | Multiskan® FC with incubator |
|--|--|------------------------------|
| Wavelength range (nm) | 340 - 850 | |
| Light source | Quartz halogen lamp 6 V/10 W | |
| Filter type | 8 position filter wheel | |
| Filters | 340 nm, 405* nm, 414 nm, 450* nm, 492 nm, 540 nm 620* nm, 690 nm filters (*standard filters) | |
| Half bandwidth of filters (nm) | 3 - 9 | |
| Linearity (96 well plate) with fast mode | 0 - 3 Abs, ±2% | |
| Linearity (96 well plate) with normal mode | 0 - 4 Abs, ±2% | |
| Read-out range (Abs) | 0 - 6 | |
| Resolution (Abs) | 0,001 | |
| Accuracy (405 nm) | ±1% (0,3 - 3 Abs) ±2% (3 - 4 Abs) | |
| Precision (405 nm) | CV = 0,2% (0,3 - 3 Abs); CV = 1,0% (3 - 4 Abs) normal mode | |
| Shaking | Linear shaking, 3 speeds | |
| Measurement time (96 well plate) | <6 s | |
| Measurement time (384 well plate) | - | <12 s |
| Display | High contrast colour display (480 x 272 dots) | |
| Software | Internal software or PC control with SkanIt® software | |
| Interfaces | USB for computer connection, USB memory stick for data export, USB for external printer (HP PCL5 type) | |
| Optional incubator | - | Included |
| Temperature range (°C) | - | Ambient +4 to 50 |
| W x D x H (mm) | 290 x 400 x 210 | |

Photometers for life science applications

Microplate reader Multiskan® GO



The Multiskan® GO is a compact, reliable microplate reader with an option available for an additional cuvette port.

The instrument can be operated as a stand-alone unit using straightforward internal software or for more demanding applications controlled by a remote PC using the SkanIt® software.

The design including a UV to visible wavelength range makes the microplate reader ideal for a wide range of photometric applications, such as nucleic acid and protein analysis, enzyme assays, cytotoxicity and cell proliferation assays as well as apoptosis assays. The option of cuvette reading port allows for optimisation of enzyme kinetics in a standard cell or TrayCell prior to transferring the kinetic protocols to the microplate format - all in the same instrument.

The SkanIt® software package supplied with the instruments provides easy assay optimisation, flexible data handling and convenient report formatting. The software also has a special remote control interface which enables easy integration with robotics and HIS/LIMS systems.

Ordering information

| Description | Cat. No. |
|------------------------------|----------|
| Multiskan GO without cuvette | 736-0559 |
| Multiskan GO with cuvette | 736-0560 |

Supplied with SkanIt® software.

- Freely selectable wavelength selection by monochromator for the demands of various assays
- Fast plate measurement and a full spectrum scan in less than 10 seconds
- Both microplate and standard 10 x10 mm cuvettes reading capability available
- Shaking and incubation up to 45 °C for temperature critical applications
- Robotic compatibility for high throughput environments

Technical specifications

| Model | Multiskan® FC | Multiskan® FC with incubator |
|--|--|---|
| Wavelength selection | Monochromator | |
| Light source | Xenon flash lamp | |
| Wavelength range (nm) | 200 -1000 with 1 nm steps | |
| Read-out range (Abs) | Up to 4 | |
| Bandwidth (nm) | 2 | |
| | Plate | Cuvette |
| Linearity (450 nm) | 0 - 2,5 Abs, ±2% (96 well plate) | 0 - 2,5 Abs, ±2% |
| Accuracy (450 nm) | 1,0% + 0,003 Abs (0 - 2,0 Abs) 2,0% (2,0 - 2,5 Abs) | 1,0% + 0,003 Abs (0 - 2,0 Abs) 2,0% (2,0 - 2,5 Abs) |
| Precision (450 nm) | CV <1,0% | CV <1,0% |
| Plate/Cuvette types | 96 and 384 well plates | Standard, micro and ultra-micro cuvettes, TrayCell |
| Measurement speed (from A1 back to A1) | 6 s with 96 well plate 10 s with 384 well plate | - |
| Shaking | Linear | - |
| Spectral scanning speed | 10 s from 200 to 1000 nm with 1 nm steps | |
| Incubation (° C) | +4 to +45 | |
| Internal user interface | 4,5" colour display USB memory stick position for data export USB port for external printer (HP PCL5 compatible) | |
| PC control | SkanIt® software using USB connection | |
| W x D x H (mm) | 285 x 430 x 260 | |
| Weight (kg) | 10,8 | |

Photometers for environmental and food industry applications

Colorimeter Spectroquant® Picco



The Spectroquant® Picco from Merck is a pocket colorimeter for water analysis.

The range covers single and multi-parameter devices that are suitable for drinking water, wastewater and swimming pool analysis.

Spectroquant® Picco colorimeters are easy to handle and are already pre-programmed for use with specific Spectroquant® test kits.

Users can choose from a series of convenient, pre-packaged cell tests and inexpensive reagent tests. Each work step is described clearly and can also be used easily by inexperienced staff.

- Simple handling
- Compact, handy and portable
- Pre-programmed methods
- Uses Spectroquant® test kits,

Items supplied:

Picco colorimeter, case, adapter for cell test kits (16 mm), 3 empty cuvettes for reagent test kits (24 mm) (not with Picco COD), 9 V battery and user manual

See the FEA catalogue W285102



Accessories for Spectroquant® Picco and Multy colorimeters

| Description | Pk | Cat. No. |
|--|----|--------------|
| Empty 24 mm cuvettes, round, with screw caps | 12 | 1.73650.0001 |
| Empty 16 mm cuvettes, round, with screw caps | 25 | 1.14724.0001 |

Technical specifications for all Spectroquant® Picco colorimeters

| | |
|--------------------------|--|
| Housing | ABS |
| Keyboard | 3 keys, polycarbonate membrane, acid and solvent-resistant, splash-proof |
| Sample chamber | Waterproof, optics temperature compensated LED |
| Power requirements | 9 V battery for approx. 40 h of use (corresponds to approx. 600 measurement series of around 4 min each) |
| Auto OFF | Automatic shut-off |
| Ambient temperature (°C) | 0 to +40 |
| Permitted rH | 30 - 90%, not condensed |
| LxWxH (mm) | 190 x 110 x 55 (without adapter), 270 x 225 x 80 (in case) |
| Weight (kg) | Device 0,4 |
| CE compliance | DIN EN 50 081-1, VDE 0839 Part 81-1 1993-03 DIN EN 50 082-2, VDE 0839 Part 82-2 1996-02 |
| Warranty | 2 years |

Ordering information

| Description | To determine | Measurement wavelength (nm) | Cat. No. |
|---|--|-------------------------------------|--------------|
| Spectroquant® Picco Cl ₂ /O ₃ /ClO ₂ /CyA/pH | Free and total chlorine, ozone, chlorine dioxide, cyanuric acid and pH | 528 (LED plus filter) | 1.73607.0001 |
| Spectroquant® Picco COD | COD | 605 (LED) and 430 (LED plus filter) | 1.73608.0001 |
| Spectroquant® Picco F | Fluoride | 620 (LED) | 1.73606.0001 |
| Spectroquant® Picco NH ₄ -N | Ammonium | 660 (LED plus filter) | 1.73602.0001 |
| Spectroquant® Picco NO ₃ -N | Nitrate | 370 (LED plus filter) | 1.73603.0001 |
| Spectroquant® Picco N (total) | Nitrogen | 370 (LED plus filter) | 1.73604.0001 |
| Spectroquant® Picco PO - P (o-phosphate and P-total) | o-Phosphate and phosphorous (total) | 660 (LED plus filter) | 1.73605.0001 |

Photometers for environmental and food industry applications

Colorimeter Spectroquant® Multy

The new Spectroquant® Multy colorimeter from Merck makes it easy to undertake photometric water analysis.

This device is pre-programmed to handle over 120 Spectroquant® reagent and cell tests and covers all major parameters of sewage water and drinking water analysis, such as COD, total nitrogen, total phosphorus and many more. You can also programme your own methods.

Even trace quantities of cyanide, arsenic, cadmium and other substances can be determined very accurately. The wide range of easily manageable Spectroquant® test kits offer an optimal solution for all users.

The device can be used with battery or mains power, so the colorimeter can be used in the field as well as in the laboratory. Method updates can be downloaded free of charge from the Internet.



- Simple to operate
- All important water, wastewater and drinking water parameters can be determined
- Portable due to the combined battery and mains connection
- Simple and free of charge updating of methods via the Internet

Items supplied:

Multy colorimeter, case, 7 rechargeable batteries, lithium battery (to ensure data preservation), battery charger, connection cable for PC, adapter and lid for 16 mm cuvettes, 3 x 16 mm round cuvettes, 3 x 24 mm round cuvettes.

Technical specifications

| | |
|--------------------------|--|
| Optics | 6 temperature compensated LEDs with interference filters, internal reference channel (dual beam technology), filters (nm) 430, 530, 560, 580, 610, 660 |
| Display | Large format graphic display |
| Output | RS232 for printer or PC connection |
| Keyboard | Acid and solvent-resistant tactile membrane keyboard with audible signal |
| Environmental conditions | Up to max. 90% relative humidity (non condensing), approx. 0 - 50 |
| Self-diagnostic | Auto check |
| Memory | Approx. 1000 data sets with date, time and registration number |
| Power | 7 NiCd rechargeable batteries (AA/Mignon), can be charged simultaneously in the device during mains operation, with overload protection |
| W x D x H (mm) | Device 265 x 195 x 70, case 440 x 370 x 140 |
| CE marked | Yes |
| Warranty | 2 years |

Ordering information

| Description | Pk | Cat. No. |
|--|----|--------------|
| Colorimeter Spectroquant® Multy | 1 | 1.73630.0001 |
| Empty 24 mm cuvettes, round, with screw caps | 12 | 1.73650.0001 |
| Empty 16 mm cuvettes, round, with screw caps | 25 | 1.14724.0001 |

Photometers for environmental and food industry applications

Photometers Spectroquant® NOVA 30A / NOVA 60 / NOVA 60A

Spectroquant® NOVA 30A

The NOVA 30A basic model from Merck is capable of running all important Spectroquant® cell tests for wastewater analysis.

The compact and mobile device can be operated with a rechargeable battery as well as mains power.

Spectroquant® NOVA 60/60A

Can be used to measure ready-to-use Spectroquant® cell tests as well as simple and affordable reagent tests. The NOVA 60 unit is a routine device with more than 170 pre-programmed test kits for Spectroquant® cell tests and 50 freely programmable methods. The NOVA 60A can also be used as a mobile analysis station. The Multi/ACHAT software makes it easier to transmit results to the PC and to programme your own methods.

- Barcode reading device for all Spectroquant® tests from Merck (automatic selection of method and then presentation of results)
- Integrated quality control - 3 modes; checking the instrument, the method and the sample matrix
- Turbidity correction
- RS232C serial interface
- For some critical tests, measurement of ppb concentration, e.g. cyanide
- Can also be used as a mobile analysis station



Technical specifications

| Model | NOVA 30A | NOVA 60A/60 |
|--------------------------|---|---|
| Kits | Most Spectroquant® tests | All Spectroquant® cell and reagent tests |
| Graphic display | 128 x 64 pixels | |
| Measurement method | 6 filters in array with the reference beam | 12 filters in array with the reference beam |
| Wavelength (nm) | 340, 445, 525, 550, 605, 690 ±2 | 340, 410, 445, 500, 525, 550, 565, 605, 620, 665, 690, 820 ±2 |
| Types of determination | Absorbance, concentration, transmission | |
| Absorbance range (E) | -0,300 to +3,200 | |
| Lamp | Tungsten halogen lamp, pre-set, no warm up time | |
| Cell compartment | 16 mm Ø cells | 16 mm Ø cuvettes and 10, 20 and 50 mm cuvettes |
| Storage capacity | 500 measured values | 1000 measured values |
| Special method functions | More than 60 pre-programmed methods | 50 freely programmable methods and more than 170 pre-programmed methods |
| Power | 110 - 130 V, 60 Hz/210 - 250 V, 50 Hz | |
| Weight (kg) | 2,8 (including battery) | |

Ordering information

| Description | Cat. No. |
|--|--------------|
| Spectroquant® NOVA 60 photometer | 1.09751.0001 |
| Spectroquant® NOVA 60A photometer | 1.09752.0001 |
| Spectroquant® NOVA 30A photometer | 1.09748.0001 |
| Multi/ACHAT II for Windows (German & English version on CD-ROM) PC software for data transfer from the Spectroquant® NOVA 30, NOVA 60 and NOVA 400 photometers. Additional control option from pH/O ₂ /conductivity meters from WTW | 1.14964.0001 |
| PC cable for Spectroquant® NOVA 30, NOVA 60, NOVA 400 photometers (for serial interface) | 1.14667.0001 |
| Halogen lamp for Spectroquant® NOVA 30 and NOVA 60 photometers | 1.09749.0001 |
| Printer cable (for serial interface) | 1.09759.0001 |
| Transport case for Nova 30 and 60 | 1.09769.0001 |

Photometers for environmental and food industry applications

Spectrophotometers Spectroquant® Pharo 100 / Pharo 300



Spectroquant® Pharo - This range of spectrophotometers combines the advantages of a system photometer - e.g. test kits ideally matched to the device - with the advantages of a spectrophotometer.

Whether you want to programme your own methods, record spectral data and kinetics or multi-wavelength measurements - all the options are available here.

These spectrophotometers have a bandwidth of 4 nm and feature excellent precision and reproducibility. Pharo spectrophotometers include barcode recognition for all Spectroquant® test kits and an automatic recognition of round or rectangular cuvettes, with no adapter required.

The large format, display, alphanumeric keyboard and user-friendly menu navigation with methods that can be updated via the Internet and on-site using a PC. This makes this spectrophotometer a versatile and flexible tool in the laboratory.

- Spectroquant® Pharo is compatible with all Spectroquant® test kits with automatic barcode recognition
- User-defined functions include recording spectral data and kinetic profiles, as well as multi-wave measurements
- Comprehensive device supported analytical quality assurance
- USB and RS232 interfaces enable the simple transfer of data to printers, PCs or in wireless mode to a USB stick
- Pharo spectrophotometers support GLP compliant work
- Flexible cell formats: Round and rectangular cuvettes (10 - 50 mm path length)

Ordering information

| Description | Cat. No. |
|--|--------------|
| Spectrophotometer Spectroquant® Pharo 100 | 1.00706.0001 |
| Spectrophotometer Spectroquant® Pharo 300 | 1.00707.0001 |
| Halogen lamp module for Spectroquant® Pharo 100 spectrophotometer | 1.00660.0001 |
| Case for Spectroquant® Pharo 100 and Pharo 300 spectrophotometers | 1.00670.0001 |
| 12 V cable for Spectroquant® Pharo 100 and Pharo 300 spectrophotometers (auto, Power Pack) | 1.00786.0001 |

* We recommend using a Power Pack rechargeable battery with a 12 V outlet for the power supply during mobile use of the Spectroquant® Pharo spectrophotometer.

Technical specifications

| Model | Pharo 100 | Pharo 300 |
|----------------------------------|--|-------------|
| Light source | Tungsten halogen lamp | Xenon flash |
| Optical system | Stabilised single beam technology | |
| Wavelength range (nm) | 320 - 1100 | 190 - 1100 |
| Measuring modes | Concentration, absorption, transmission, multi-wavelengths, scans + kinetics | |
| Spectral bandwidth (nm) | 4 | |
| Wavelength accuracy (nm) | ±1 | |
| Absorbance accuracy | 0,003 E at <0,600 E 0,5% of the measured value for 0,600 <A >2,000 | |
| Photometric range (E) | ±3,3 | |
| Scan | 1 nm steps with a selectable wavelength range | |
| Cuvettes | 16 mm round, 10/20/50 mm rectangular with automatic recognition | |
| Interface | RS232, USB-A, USB-B | |
| Data storage | 1000 individual measured values, 4 MB for scans and kinetics | |
| Protection class | IP31 and drain in the cell compartment | |
| Power | 100 - 240 V/50 - 60 Hz | |
| W x D x H (mm) | 404 x 314 x 197 | |
| Weight (kg) without power supply | 3,7 | |

Photometers for environmental and food industry applications

Spectrophotometer AquaMate™ plus for water analysis



A powerful single beam spectrophotometer from Thermo Scientific with quartz coated optics for water analysis.

The methods and results are shown on a VGA quality LCD display.

The USB interface allows a wide range of methods and data to be stored. The device has an internal memory that can store the 20 most frequent water analysis methods.

It is also supplied with the calibration data and factors from Merck Spectroquant®, Hach®, CHEMetrics® and Dr. Lange test kits. The Uvcalc™ internal calculation software means data can be evaluated simply and flexibly.

- Methods can be password protected
- Easy loading of standard methods and updating of existing methods
- Sipper systems and numerous other accessories available as options

Ordering information

| Description | Cat. No. |
|--|----------|
| Spectrophotometer AquaMate™ plus VIS, 325 - 1100 nm, EU and UK plug | 705-0951 |
| Spectrophotometer AquaMate™ plus UV/VIS, 190 - 1100 nm, EU and UK plug | 705-0950 |

Software details available on Page 23

Accessories

| Description | Cat. No. |
|--|----------|
| 7x cell changer, automatic | 635-2165 |
| Test tube/cell holder combination | 634-1028 |
| Variable long cell holder for rectangular cuvettes 1 - 50 mm layer thickness | 705-0608 |
| Long cell holder for 100 mm rectangular cuvettes | 705-0329 |
| Long cell holder for 100 mm cylindrical cuvettes | 705-0327 |
| Cell holder for Hach® square 1" cells and AccuVac® ampoules | 705-0330 |
| Thermostatable single cell holder, 1 to 50 mm layer thickness | 634-1027 |
| Single cell holder, Peltier tempered | 634-1030 |
| SuperSipper | 634-1029 |
| MiniSipper | 635-2163 |

Technical specifications

| Model | AquaMate™ plus VIS | AquaMate™ plus UV/VIS |
|-----------------------------------|--|------------------------|
| Light source | Tungsten | Deuterium and tungsten |
| Detector | Silicon photodiode | |
| Optical system | Single beam | |
| Range (nm) | 315 - 1100 | 190 - 1100 |
| Spectral data scan speed (nm/min) | 1 - 380 | |
| Scan intervals (nm) | 0,2; 0,5; 1,0; 2,0; 4,0; 10,0 | |
| Spectral bandwidth (nm) | 2 | |
| Wavelength accuracy (nm) | ±1,0 | |
| Absorbance accuracy | 0,005 A at 1 A | |
| Photometric range | -0,1 to +200% T; -0,3 to +3,0 A, 0 - 9999 C | |
| Display | VGA LCD | |
| Keyboard | Membrane keypad | |
| Interfaces | RS232C, USB for printer and data storage | |
| Memory | Internal: Up to 30 methods and data set calibrations or external: Unlimited methods and data sets via USB on memory stick | |
| Standard cell holder supplied. | Test tube (max. Ø 18 mm)/rectangular cell holder combination for tubes with max. Ø 18 mm, total length of 125 mm and rectangular cuvettes with 10 mm path length. Variable rectangular cell holder for cuvettes from 1 - 50 mm path length. Cell holder for Hach®, 1" rectangular cuvettes and AccuVac® tubes. | |
| Software | Wavelength scan with peak search, absorbance, transmission, concentration, kinetic, multi-component analyses, multi-wavelength measurements (up to 20), standard curves (with up to 20 standards), time scan | |
| Optional software | VISION software e.g. VISION/ite™, EnzLab, ColorCalc | |
| Power (V) | 100 - 240 | |
| W x D x H (mm) | 455 x 395 x 215 | |
| Weight (kg) | 10 | |
| Warranty | 3 years | |

Universal photometers - single beam devices

Spectrophotometers VWR Collection V-1200



The V-1200 spectrophotometer is a basic VIS unit attractively priced with robust housing.

It offers all standard methods including photometry mode (absorption, transmission) and standard curves (with standards or coefficient mode). Using a standard sample solution, you can get a standard curve on the large LCD screen by the locally controlled software, and then print the curve through the parallel port.

The unit is ideal for educational establishments and simple daily quality control.

- Large LCD screen (128 x 64 dots) can display a total of 50 groups of data, 3 groups per screen. It can also display standard curve and the curve equations
- Internal memory can save 200 sets of data and 200 standard curves which is convenient for check and reload
- Pre-aligned design makes it convenient to change lamps
- Large sample compartment which can accommodate 5 - 100 mm path length cuvettes with optional holders
- Variety of optional accessories are available

Ordering information

| Description | Cat. No. |
|--|----------|
| V-1200 spectrophotometer including 4 glass cells, power cords UK/EU/SW, manual | 634-6000 |

Accessories


| Cell holders for | Cat. No. |
|-----------------------------------|----------|
| 4 rectangular cells, 10 to 50 mm | 634-6004 |
| 4 rectangular cells, 10 to 100 mm | 634-6005 |
| Cylindrical cells | 634-6006 |
| Micro cells | 634-6008 |
| Test tubes | 634-6009 |
| Solid samples | 634-6011 |

Spares

| Description | Cat. No. |
|-------------------------|----------|
| Halogen lamp, 12 V/20 W | 634-6037 |
| Thermal printer | 634-6039 |

Technical specifications

| | |
|--------------------------|--|
| Optical system | Single beam |
| Wavelength range (nm) | 325 - 1000 |
| Band width (nm) | 4 |
| Stray light | ≤0,2%T |
| Photometric range | 0 - 200% T, -0,3 to +3,0 A, 0 - 9999 C |
| Wavelength accuracy (nm) | ±2 |
| Photometric accuracy | ≤0,5%T or 0,005 A @ 1 A |
| Stability | 0,002 A/h @ 500 nm |
| Memory | 200 results & 200 standard curves |
| Languages | EN, FR, DE, SP |
| Display | 128 x 64 dots graphic LC display |
| Interface | USB, printer port |
| Methods | Photometry (absorbance, transmittance, concentration by linear function or standard curve-up to 9 standards) |
| W x D x H (mm) | 490 x 360 x 210 |
| Weight (kg) | 12 |
| Warranty | 2 years |



| Quantitation | | | |
|-----------------------------------|---------|-------|-------|
| Date & Time: mm-dd-yyyy, hh:mm:ss | | | |
| Model: V-1200 | | | |
| Serial No.: VECXXXXXXXX | | | |
| Firmware Version: 2.5.0 | | | |
| VWR International bvba. | | | |
| C=1.000*A+1.000 | | | |
| r=1.0000 | | | |
| No. | WL.(nm) | Abs | Conc. |
| 1 | 500.0 | 0.120 | 1.120 |
| 2 | 500.0 | 0.127 | 1.127 |
| 3 | 500.0 | 0.121 | 1.121 |
| End. | | | |

Universal photometers - single beam devices

Spectrophotometer GENESYS™ 20 VIS

GENESYS™ 20 from Thermo Scientific

This single beam photometer is ideal for taking routine measurements in teaching laboratories, process control and production. It is reliable, robust, accurate and can also be supplied with an integrated printer.

A wide range of simple cuvette holders, e.g. for COD vials, 50 mm cuvettes, filters and test tubes, are also available for installation in the instrument.

- Splash-proof membrane keyboard for easy cleaning
- Clear 2 line display
- Easy-to-use



Accessories

| Description | Cat. No. |
|--|----------|
| Printer paper roll for the internal printer for GENESYS™ 20 | 634-0502 |
| Standard cell holder for cuvettes with a 10 mm layer thickness | 634-0388 |
| Replacement halogen lamp for GENESYS™ 20 | 634-1000 |
| 1 inch round cuvette holder | 634-1673 |
| Standard filter set | 634-2113 |
| Single cell holder (suitable for a standard platform and 6x changer) | 634-1669 |

Other accessories are available on request.

Technical specifications

| | |
|------------------------------|---|
| Light source | Tungsten lamp |
| Optical system | Single beam |
| Measuring range (nm) | 325 - 1100 |
| Spectral bandwidth (nm) | 8 |
| Accuracy (nm) | ±2,0 |
| Photometric range | -0,1 to +2,5 A; 0 - 25% T; 0 - 1999 C, absorption, transmission, concentration, factors |
| Display | LCD with 2 lines and 20 characters |
| Standard cell holder | Holder for cuvettes with a 10 mm layer thickness and tubes |
| Test storage | 1 set of parameters in a ROM |
| Storage capacity | None |
| Printer, internal (optional) | 20 columns |
| W x D x H (mm) | 300 x 330 x 190 |
| Weight (kg) | 4,5 |
| Warranty | 3 years |
| Accessories included | Dust cover, user manual, 5 disposable cuvettes |



Ordering information

| Description | Cat. No. |
|--|----------|
| GENESYS™ 20, EU plug | 634-1021 |
| GENESYS™ 20, UK plug | 634-0385 |
| GENESYS™ 20 with integrated thermal printer, EU plug | 634-1022 |
| GENESYS™ 20 with integrated thermal printer, UK plug | 634-0387 |


See page 25 for details of the software packages

Universal photometers - single beam devices

Spectrophotometers VWR Collection UV-1600 PC



| No. | WL. | Abs. |
|-----|-------|-------|
| 1 | 230.0 | 0.001 |
| 2 | 340.0 | 0.000 |
| 3 | 450.0 | 0.002 |
| 4 | 540.0 | 0.000 |
| 5 | 620.0 | 0.003 |



UV-1600PC

The UV-1600PC is a full range UV/VIS spectrophotometer 190 to 1100 nm. Additional features include kinetics and time scan. Remote control software enables you to run wavelength scans and multi-wavelength measurements.

- Full wavelength range from 190 to 1100 nm
- Internal memory can save 200 sets of data and 200 standard curves which is convenient for check and reload
- Pre-aligned design makes it convenient to change lamps
- Large sample compartment, can accommodate 5 - 100 mm path length cuvettes and optional holders up to electronic thermostatted cell holder, sipper and more
- Variety of optional accessories available

Technical specifications

| | |
|--------------------------|---|
| Optical system | Single beam |
| Wavelength range (nm) | 190 - 1100 |
| Band width (nm) | 4 |
| Stray light | ≤0,05% T @ 220 nm & 360 nm |
| Photometric range | 0 - 200% T, -0,3 to +3,0 A, 0 - 9999 C |
| Wavelength accuracy (nm) | ±0,5 |
| Photometric accuracy | ≤±0,5% T or 0,005 A @ 1 A |
| Stability | 0,002 A/h @ 500 nm |
| Memory | 200 results & 200 standard curves |
| Languages | EN, FR, DE, SP |
| Display | 128 × 64 dot matrix LCD |
| Interface | USB, Parallel |
| Methods | Built-in: Photometry (absorbance, transmittance, concentration by linear function or standard curve-up to 9 standards), kinetics/time scan. Methods by PC remote control: Photometry (absorbance, transmittance, concentration by linear function or standard curve-up to 20 standards, multi-wavelength measurements – up to 20 wavelengths), DNA/protein analysis, kinetics/time scan, wavelength scan. |
| Power | AC 110/220 V, 50/60 Hz |
| W x D x H (mm) | 490 × 360 × 240 |
| Weight (kg) | 14 |
| Warranty | 2 years |

Ordering information

| Description | Cat. No. |
|---|----------|
| UV-1600PC spectrophotometer including 4 glass cells, 2 quartz cells, UK, EU, SW, USB cables, remote control software and manual | 634-6001 |

Universal photometers - single beam devices


Spectrophotometers VWR Collection UV-3100 PC



UV-3100PC

The UV-3100 PC represents the top end of our 3 spectrophotometers with a 2 nm spectral bandwidth for a higher resolution and peak sensitivity. A large graphic display helps you to get all results and data including standard curves at a glance.

All methods are available at built-in/stand-alone software and on the remote control software for a flexible use of the instrument.



- Internal memory can save 200 sets of data and 200 standard curves which is convenient for check and reload
- Pre-aligned design makes it convenient to change lamps
- Large sample compartment can accommodate 5 - 100 mm path length cuvettes and optional holders up to electronic thermostatted cell holder, sipper and more
- Variety of optional accessories available
- Large graphic LCD display; 320 x 240 dots

Spares for UV-1600 PC and UV-3100 PC

| Description | Cat. No. |
|-------------------------|----------|
| Halogen lamp, 12 V/20 W | 634-6037 |
| Deuterium lamp | 634-6038 |
| Thermal printer | 634-6039 |

Accessories for UV-1600 PC and UV-3100 PC

| Description | Cat. No. |
|---|----------|
| Cell holders | |
| For 4 rectangular cells, 10 to 50 mm | 634-6004 |
| For 4 rectangular cells, 10 to 100 mm | 634-6005 |
| For cylindrical cells | 634-6006 |
| Thermostatable for 1 x 10 mm rectangular cell | 634-6007 |
| For micro cells | 634-6008 |
| For test tubes | 634-6009 |
| Automatic 8 position for 10 mm rectangular cells | 634-6010 |
| For solid samples | 634-6011 |
| Thermostatable for 4 x 10 mm rectangular cells | 634-6012 |
| Peltier thermostated for 1 x 10 mm rectangular cell | 634-6034 |
| Sippers | |
| Sipper unit | 634-6035 |
| Sipper unit with temperature control | 634-6036 |

Technical specifications

| | |
|--------------------------|--|
| Optical system | Single beam, grating 1200 lines/mm |
| Wavelength range (nm) | 190 - 1100 |
| Band width (nm) | 2 |
| Stray light | ≤0,05% T @ 220 nm & 360 nm |
| Photometric range | 0 - 200% T, -0,3 to +3,0 A, 0 - 9999 C |
| Wavelength accuracy (nm) | ±0,5 |
| Photometric accuracy | ≤±0,5% T or 0,005 A @ 1 A |
| Stability | 0,002 A/h @ 500 nm |
| Scan speed (nm/min) | Max. 300 |
| Memory | 200 results & 200 standard curves |
| Languages | EN, FR, DE, SP |
| Display | 320 x 240 graphic LCD |
| Interface | USB, parallel |
| Methods | Built-in: Photometry (absorbance, transmittance, concentration by linear function or standard curve-up to 9 standards), multi-wavelength photometry, kinetics/time scan, Methods by PC remote control: Photometry (absorbance, transmittance, concentration by linear function or standard curve-up to 20 standards, multi-wavelength measurements – up to 20 wavelengths), DNA/protein analysis, kinetics/time scan, wavelength scan |
| Power supply | AC 110/220 V, 50/60 Hz |
| W x D x H (mm) | 490 x 360 x 240 |
| Weight (kg) | 14 |
| Warranty | 2 years |

Ordering information

| Description | Cat. No. |
|---|----------|
| UV-3100PC spectrophotometer including 4 glass cells, 2 quartz cells, UK, EU, SW, USB cables, remote control software and manual | 634-6002 |

Universal photometers - single beam devices

Photometer GENESYS™ 10S VIS



GENESYS™ 10S

This series of instruments offer excellent value for money. They also provide robustness with ease-of-use.

The efficient optical configuration delivers high performance with a minimum number of optical surfaces in a true monochromator design. The xenon lamp provides instant UV/VIS measurements and is guaranteed for 3 years of continuous use.

The GENESYS™ 10S VIS uses a tungsten lamp and single detector to support routine measurements in the visible range.

- Spectrum scan function across the complete measuring range from 190 - 1100 nm at 10 - 4200 nm/min scan speed
- USB printer interface for text and graphics in HP PCL format
- Automatic correction of layer thickness tolerances with up to 6 cuvettes
- SmartStart™ software offers a simple selection of standard methods when used for the first time

Ordering information

| Description | Cat. No. |
|---|----------|
| GENESYS™ 10S VIS, power leads with EU and UK plug | 634-0593 |
| GENESYS™ 10S VIS with internal printer and power leads with EU and UK plug | 634-0594 |
| GENESYS™ 10S UV/VIS, power leads with EU and UK plug | 634-0595 |
| GENESYS™ 10S UV/VIS with internal printer and power leads with EU and UK plug | 634-0596 |

Supplied with a 6 position cell changer, single cell holder, spare fuses and USB memory device.

Technical specifications

| Model | GENESYS™ 10S VIS | GENESYS™ 10S UV/VIS |
|--------------------------|--|---|
| Lamp source, lifetime | Tungsten halogen (1000 h) | Xenon flash (5 years typically) |
| Optical system | Single beam | Dual beam internal reference detector |
| Range (nm) | 325 - 1100 | 190 - 1100 |
| Spectral bandwidth (nm) | 5,0 | 1,8 |
| Wavelength accuracy (nm) | ±1,0 | |
| Scan speed (nm/min) | 10 - 4200 | |
| Photometric range | Up to 3,0 A at 340 nm | Up to 3,5 A at 260 nm |
| Photometric accuracy | 0,5% or ±0,005 A, whichever is greater up to 2 A | ±0,005 A at 1,0 A 0,010 A K2Cr2O7 |
| Stray light | <0,1%T at 340, and 400 nm | <0,08% T at 220, 340 nm (NaI, NaNO ₂) <1,0% 198 nm (KCl) |
| Display | Graphic 320 x 240 pixel backlit LCD (91 x 71 mm) | |
| Standard cell holder | Integral 6 position cell changer, single cell holder | |
| Keyboard | Membrane keypad | |
| Printer | Optional, 40 column, internal, graphic | |
| Interfaces | USB type A port for USB memory device (front panel) USB type B port for optional PC (rear panel) USB type A port for external printer (rear panel) | |
| Power (V) | 100 - 240, selected automatically | |
| W x D x H (mm) | 300 x 400 x 250 | |
| Weight (kg) | 8,6 | |

Universal photometers - single beam devices

GENESYS™ 20/10S/10S BIO/Aquamate plus photometers software

Software for GENESYS™ 20/10S/10S BIO/AquaMate™ plus photometers

It contains all common measuring modes; standard photometry (Fixed), spectrum scan (Scan), determination of concentration (Quant) and kinetic (Rate).

Ordering information


| Description | Cat. No. |
|---|----------|
| VISION/ite™ software, CD-ROM and user manual | 634-1017 |
| PC connection cable (DB-9) for GENESYS™ 6/10/20 | 634-1016 |

- Fixed mode carries out single and multi-wavelength measurements (up to 31), for rapid identification/ changing of concentrations
- Scan mode records a partial or full spectrum for your sample
- Rate mode measures your sample at a selected wavelength depending on the time
- Quant mode determines the sample concentration using either a standard curve or a factor
- Other options include measuring up to 100 samples; up to 13 graphics can be overlaid
- Peak identification, reference wavelength, 4 curve models, curve with up to 20 standards, sample data reporting with pre-programmed units, automatic data logging and export

Photometers Ultrospec 5300/6300 pro

High performance instruments featuring 1 nm bandwidth for precision measurements in analytical and research laboratories.

An 8 position sample changer comes as standard. Accessories such as a sipper option are also available.



Ordering information

| Description | Colour | Cat. No. |
|---|---------|------------|
| Ultrospec 5300 PC-controlled instrument with SWIFT II software supplied | Gold | 80-2117-56 |
| | Classic | 80-2117-70 |
| | Yellow | 80-2117-71 |
| | Plum | 80-2117-72 |
| | Apple | 80-2117-73 |
| Ultrospec 6300 stand-alone instrument with built-in methods | Gold | 80-2117-55 |
| | Classic | 80-2117-60 |
| | Yellow | 80-2117-61 |
| | Plum | 80-2117-62 |
| | Apple | 80-2117-63 |

- Built-in applications for nucleic acid and protein quantitation, cDNA software (measures Cy™ 3, Cy™ 5 and fluorescein) for microarray probes and PCR products; enzyme kinetics, wavelength scanning, substrate quantitation, and standard curves
- Wide wavelength range 190 to 1100 nm in 0,1 nm data steps for flexible application development
- Wavelength accuracy $\pm 0,5$ nm
- Absorbance range: -3,000 to +3,000 A, $\pm 0,5\%$ or $\pm 0,003$ A (at 546 nm), whichever is greater
- Reference beam compensation
- Rapid scan feature
- GLP self-diagnostics
- IQ/OQ certification test plans
- Compatible with 21 CFR Part 11
- Pharmacopoeia compliant with qualification and performance verification logbook as standard

Universal photometers - single beam devices

UV/VIS spectrophotometer HITACHI U-5100



Ordering information

| Description | Cat. No. |
|---|----------|
| U-5100 UV/VIS spectrophotometer | 634-0757 |
| Autosipper | 634-0758 |
| Single cell holder for 10 mm rectangular cells | 634-0759 |
| Cell holder for rectangular cells, 10 to 100 mm | 634-0760 |



HITACHI's U-5100 spectrophotometer has a small footprint and the functionality to accommodate a wide range of applications.

Designed with the environment in mind, the U-5100 has a long life xenon flash lamp which eliminates the need for periodic lamp replacement. The utilisation of a xenon flash lamp also reduces power consumption and noise - once the unit is switched on, sound is minimal during measurements and silent in standby mode.

The user interface has been designed to input and display clear information with a large 6 inch LCD mounted display at an easy-to-read angle above a soft keypad. The keys used for character input, necessary for assigning a file name, are designed to be as easy-to-use as texting using a mobile phone.

The standard automatic 6 cell changer makes daily work more efficient. Standard photometry, time scan, wavelength scan, kinetics and performance validation functionality are all included as standard in the menu driven internal software.

- Low noise xenon flash lamp uses less power
- Soft keypad
- Large 15 cm display for easy reading and programming
- PC control with UV Solutions® software (optional)
- Automatic 6 position cell changer as standard

Technical specifications

| | |
|---------------------------------|---|
| Optics | Seya-Namioka mount monochromator, ratio beam |
| Detector | Silicon photodiode |
| Light source | Silicon flash lamp |
| Spectral band pass (nm) | 5 |
| Wavelength range (nm) | 190 – 1100 |
| Wavelength accuracy (nm) | ±1 |
| Photometric range | –3,000 to +3,000 A/0 – 300%T/0,000 – 9999 conc. |
| Photometric accuracy | ±0,003 A (0 to 0,5 A)/±0,005 A (0 to 1,0 A) |
| Photometric repeatability (Abs) | ±0,002 A (0 to 1,0 A) |
| Stray light | ≤0,07% at 220 and 340 nm |
| Baseline stability | 0,0007 A/h (260 nm) |
| Noise level | 0,0002 A or less (RMS (260 nm), 0 A) |
| Baseline flatness | ±0,010 A (200 – 950 nm) |
| Wavelength scan speed (nm/min) | 40, 100, 200, 400, 800, 1200, 2400 |
| Display | 15 cm LCD with backlight |
| Languages | EN, DE, SP |
| Performance test | Validation function |
| Standard cell holder | Automatic 6 cell turret (cuvettes not supplied) |
| Interfaces | Centronics interface (parallel port connection); USB for PC control |
| Internal memory | 50 methods, 30 sets of data |
| Power | 100 – 240 V; 50/60 Hz; 60 W |
| W x D x H (mm) | 355 x 425 x 235 |
| Weight (kg) | 13 |

Universal photometers - single beam devices

UV/VIS spectrophotometers UviLine 9100 / 9400

UviLine 9100/9400 from SI Analytics

A high performance, single beam spectrophotometer, with choice of VIS (9100) or UV/VIS (9400) wavelength ranges.

In addition to the standard spectrometer functions of absorbance and transmission, the instrument has a wide range of menu driven evaluations and therefore has applications across many disciplines. Furthermore the user can easily programme and store more than 100 of their own methods. A large graphic display provides online spectra and together with the simple to use keyboard guides the user with menu navigation.

- Complete functionality for scanning with online graphics, kinetics and multi-wavelength analytics
- USB Master interface (USB-A) enables connection with USB memory sticks, printers with USB interface and external PC keyboards
- Automatic wavelength calibration and compensation of ambient light
- Extensive range of evaluation functions such as min./max. recognition, add and subtract spectrum etc.
- GLP compliant with user login and 3 levels of operation



Ordering information

| Description | Cat. No. |
|--------------|----------|
| UviLine 9100 | 634-5000 |
| UviLine 9400 | 634-5001 |

Accessories

| Description | Cat. No. |
|--|----------|
| Replacement halogen lamp for UviLine 9100 | 634-0597 |
| 5+1 automated cell changer | 634-5002 |
| Cell holder, single, 10 mm | 634-0598 |
| Cell holder, single, thermostatable, 10 mm | 634-5004 |
| Sipper SZ2150 | 634-5005 |
| Cell holder, single, thermostatable (Peltier), 10 mm | 634-5006 |

Technical specifications

| Model | UviLine 9100 | UviLine 9400 |
|-------------------------|---|--------------|
| Accuracy (nm) | ±1 | |
| Display | Graphic backlit display, 320 x 240 pixels | |
| Interfaces | 1 x USB-A; 1 x USB-B; 1 x RS232C | |
| Keypad | Alphanumeric, function and cursor direction soft keys | |
| Light source | Tungsten-Halogen | Xenon |
| Method storage | >100 (extendable via USB) | |
| Optical system | Single beam | |
| Photometric accuracy | 0,3% or ±0,003 A (from 0 - 0,6 A) | |
| Photometric range | -3,3 to +3,3 A | |
| Photometric stability | <1% at 2 A between 340 - 900 nm | |
| Power | 110 - 220 V, 50/60 Hz | |
| Range (nm) | 320 - 1100 | 190 - 1100 |
| Spectral bandwidth (nm) | 4 | |
| Stray light | <0,1% at 220, 340 and 400 nm | |
| Weight (kg) | 4 | |
| W x D x H (mm) | 404 x 314 x 197 | |

Supplied with a thermostatable single cell holder (10 mm path length).

Universal photometers - single beam devices

Universal photometers - single beam devices Jenway 7300 series



The 73 series spectrophotometers are designed to provide the user with a compact, easy to operate instrument for use in education and general QC laboratories.

- Icon driven software for easy and intuitive navigation
- Small footprint with graphic display in the compartment lid
- Autologging capabilities and optional internal printer
- Press-to-read xenon lamp (7305 and 7315) to extend xenon lamp life
- Extensive range of accessories available

Models 7310 and 7315

These models have more advanced features including:

- Scanning capability
- Kinetic measurements with real-time graphic display
- Supervisor security functions to protect methods/options and operator ID
- Quantitation, using up to 6 standards with curve fit
- Results and method saving to USB memory stick

Ordering information

| Description | Cat. No. |
|--------------------------------------|----------|
| 7300 spectrophotometer, VIS range | 664-0067 |
| 7305 Spectrophotometer, UV/VIS range | 664-0066 |
| 7310 spectrophotometer, VIS range | 664-0065 |
| 7315 Spectrophotometer, UV/VIS range | 664-0064 |

Accessories

| Description | Cat. No. |
|---|----------|
| Sipper pump | 664-0063 |
| Peltier pump | 664-0062 |
| Sipper/Peltier pump | 664-0060 |
| 8 cell automatic turret | 664-0061 |
| 10 x 10 mm path length cell holder | 664-0068 |
| 16/24 mm test tube holder | 664-0059 |
| 10 to 100 mm path length cell holder | 664-0058 |
| Micro cell holder with reduced aperture | 664-0057 |
| Single cell holder, water heated | 664-0056 |
| Lamp, tungsten halogen for 6300/7300 | 634-0083 |
| Internal printer for 7300 series | 634-3009 |

Technical specifications

| Model | 7300 | 7305 | 7310 | 7315 |
|-------------------------------|---|------------|---|------------|
| Light source | Tungsten halogen | Xenon | Tungsten halogen | Xenon |
| Optical system | Single beam with silicon photodiode | | | |
| Wavelength range (nm) | 320 - 1000 | 198 - 1000 | 320 - 1000 | 198 - 1000 |
| Spectral bandwidth (nm) | 5 | | | |
| Wavelength accuracy (nm) | ±2,0 nm | | | |
| Wavelength resolution (nm) | 1 | | | |
| Wavelength repeatability (nm) | ±0,5 nm | | | |
| Photometric range | 0 - 199,9% T; -0,300 to +2,500 A | | | |
| Photometric accuracy | ±1,0% T; ±0,01 A at 1,000 A | | | |
| Photometric resolution | 0,1% T; 0,001 A | | | |
| Measuring modes | Absorbance, % transmittance and concentration | | Absorbance, % transmittance, concentration; scanning, kinetics and quantitation | |
| Calibration | Blank with a single standard or factor | | | |
| Display | Graphic LCD | | | |
| Interfaces | RS232 | | RS232 + USB port | |
| W x D x H (mm) | 275 x 400 x 220 | | | |
| Weight (kg) | 6 | | | |

Delivery information

Spectrophotometers fitted with 10 x10 mm cuvette holder, supplied with a pack of 100 disposable cuvettes, instruction manual, universal power supply, PC software on CD ROM and interface cable. Models 7310 and 7315 are supplied with USB memory stick.

Universal photometers - double beam devices

Split beam spectrophotometer 6700 series



A range of scanning spectrophotometers with a wide range of applications that have been designed for environments where speed, sample flexibility and regulatory conformance are required.

Scanning at an impressive 1500 nm/min even with data collection at 0,1 nm intervals. This enables zooming to higher resolutions without the need to re-scan. The large colour touch screen display offers simple, intuitive operation and clear representation of results. The sealed QWheel™ allows for fine cursor control or scrolling. The range features significant data storage capability with the option of a removable SD and SD/USB memory cards to transfer data to another instrument or a PC. Transfer of data is also possible via the included software.

Ordering information

| Description | Cat. No. |
|---|----------|
| 6700 VIS with automatic 8 cell changer 634-3001 | 634-3001 |
| 6700 VIS with automatic 8 cell changer and internal printer | 634-0583 |
| 6700 VIS with single cell holder 634-0583 | 634-0533 |
| 6700 VIS with single cell holder and internal printer | 634-3003 |
| 6705 VIS with automatic 8 cell changer 634-3003 | 634-0584 |
| 6705 VIS with automatic 8 cell changer and internal printer | 634-0534 |
| 6705 VIS with single cell holder 634-0584 | 635-3005 |
| 6705 VIS with single cell holder and internal printer | 634-0585 |
| 6715 VIS with automatic 8 cell changer 634-3005 | 634-0535 |
| 6715 VIS with automatic 8 cell changer and internal printer | 634-3001 |
| 6715 VIS with single cell holder 634-0585 | 634-0583 |
| 6715 VIS with single cell holder and internal printer | 634-3003 |

Accessories

| Description | Cat. No. |
|---------------------------------------|----------|
| Sipper system for 6700 series | 634-3006 |
| Peltier for 6700 series | 634-3007 |
| Water jacketed cell holder | 634-3008 |
| Internal printer for 6700 series | 634-3009 |
| Sipper/Peltier system for 6700 series | 634-3010 |

- Colour LCD with touch screen interface
- Individual operator passwords and different levels of access provide secure multi-user operation
- Kinetics for timed measurement and enzyme reaction monitoring
- Plug-in accessory modules for a range of sample options
- Software supports compliance with CFR21 Part 11

Technical specifications

| Model | 6700 | 6705 | 6715 |
|-------------------------------------|--|---------------|-------|
| Light source | Tungsten | Halogen | Xenon |
| Optical system | Sealed, MgF coated, split beam | | |
| Spectral bandwidth (nm) | 4 | | 1,5 |
| Stray light (nm) | <0,1% at 340 | <0,05% at 220 | |
| Wavelength range (nm) | 320 to 1100 | 190 to 1100 | |
| Wavelength resolution/accuracy (nm) | 0,1/±1,0 | | |
| Wavelength repeatability (nm) | ±0,2 | | |
| Photometric range | -0,300 to +3,000 A & 0 to 199,9% T | | |
| Photometric resolution | 0,001 A & 0,1% T | | |
| Photometric accuracy (A) | ±0,005 at 1 | | |
| Photometric stability | <0,001 A per h | | |
| Quantitation range | -99 999 to +99 999 | | |
| No. of standards | 20 with up to 5 replicates of each | | |
| Curve fit algorithms | Linear, quadratic and cubic functions | | |
| Multi-wavelength data points | Up to 4 wavelengths | | |
| Calculations | Difference and ratio | | |
| Kinetics time limits (s) | 0 to 9999 | | |
| Kinetics calibration | Standard or factor | | |
| Scan speed | 1500 nm/min at 0,1 nm data steps | | |
| Post scan analysis | Peak/valley pick, peak ratios, area, zoom, wavelength table, derivatives, smoothing, overlay | | |
| Configuration | Secure multi-user or free access for up to 10 users | | |
| No. of methods | >1000 on internal flash memory or removable media | | |
| Results storage | >1000 on internal flash memory or removable media | | |
| Removable media | MM/SD memory card or SD/USB memory card | | |
| Interfaces | USB, Centronics, analogue | | |
| PC software | Supplied on CD ROM with USB interface cable | | |
| Power | 100 to 230 V, AC 50/60 Hz (UK, Euro or US power leads supplied) | | |
| WxHxD (mm) | 490 x 390 x 220 | | |
| Weight (kg) | 7,5 | | |

Universal photometers - double beam devices

UV/VIS spectrophotometer U-2900 / U-2910



Reliable and accurate spectrophotometers from HITACHI with dual beam technology, for use in quality control and research for many different applications including life science and development of new materials.

The U-2900 has been designed for direct use without a PC as a stand-alone device. The keyboard and the extraordinary large graphic display make it easy to enter method data and convenient to display results. Methods and results can also be saved on a USB memory stick in the device format or as a text/table file. The device also has the option of being controlled externally via a PC using UV Solutions® software. The U-2910 was developed exclusively for connection to a PC and has therefore dispensed with the keyboard and display, although it has the same equipment otherwise.

Software

| Description | Cat. No. |
|--|----------|
| UV Solutions® software | 634-0724 |
| Optional package for colour measurements | 634-0725 |
| Report generator | 634-0726 |
| Nucleic acid measurement | 634-0728 |
| GLP/GMP programme | 634-0729 |

- Complies with the specifications of the European Pharmacopoeia Ph. Eur.
- Large, adjustable LCD display (10,4") and user-friendly keyboard (U-2900)
- Full scans across the complete wavelength range of 190 to 1100 nm with a selectable scan speed of 10 to 3600 nm/min
- Validation and self-test functions
- USB interface for storing methods and data (on a memory stick)

Accessories

| Description | Cat. No. |
|---|----------|
| Autosipper | 634-0720 |
| Autosipper with Peltier temperature control from 20 to 40 °C, for 10 mm cells | 634-0721 |
| AS-1010 autosampler (autosipper required) | 634-0723 |
| Automatic 6 cell holder with Peltier temperature control, 20 to 40 °C | 634-0722 |
| Single cell holder with Peltier temperature control, 10 to 60 °C and magnetic stirrer for 10 mm cells | 634-0050 |
| Thermostatable single cell holder, thermostat required | 634-0059 |
| Tandem cell holder for max. 3 cells on sample and reference position | 634-0062 |
| Single rectangular long cell holder for 10, 20, 30, 40, 50, 100 mm | 634-0055 |
| 5 cell changer, manual | 634-0058 |
| 4 position rectangular long path cell holder for 10, 20, 30, 40, 50, 100 mm | 634-0054 |
| Glass filter holder | 634-0057 |
| Single long cell holder for cylindrical cells | 634-0056 |
| Film holder | 634-0060 |

Technical specifications

| Model | U-2900 | U-2910 |
|--------------------------|--|-----------------|
| Light source | Tungsten, deuterium lamp | |
| Optical system | Dual beam | |
| Measuring range (nm) | 190 - 1100 | |
| Scan speed (nm/min) | 10, 100, 200, 400, 800, 1200, 2400, 3600 | |
| Spectral bandwidth (nm) | 1,5 | |
| Wavelength accuracy (nm) | ±0,3 (at 656,1; 486,0) | |
| Photometric range | 0 - 300% T; -3 to +340 A | |
| Display | Colour LCD, backlit | - |
| Keyboard | Robust laboratory keyboard | - |
| Test storage | Internal memory or USB memory stick | - |
| Measurement methods | Absorption, transmission, ratio, concentration, wavelength scan, time scan, multi-wavelength measurement ratio calculations from A260/A280, validation functions and automatic calibration and self-test functions | |
| Interfaces | RS232C for PC, Centronics for printer, USB port for memory stick | |
| Power (V) | 100 - 240 | |
| W x D x H (mm) | 500 x 605 x 283 (folded down LCD) | 500 x 605 x 241 |
| Weight (kg) | 31 | 29 |

Ordering information

| Description | Cat. No. |
|--|----------|
| U-2900 spectrophotometer | 634-0718 |
| U-2910 spectrophotometer including UV Solutions® software (PC required) only EU plug available | 634-0719 |

Universal photometers - double beam devices

Double beam UV/VIS spectrophotometers U-3900 / U-3900H



Accessories

| Description | Cat. No. |
|--|----------|
| Autosipper | 634-0577 |
| Autosipper with Peltier temperature control | 634-0578 |
| AS-1010 autosampler, autosipper required | 634-0723 |
| Single cell holder with Peltier temperature control, 10 to 60 °C, for 10 mm cells | 634-0050 |
| Single cell holder with Peltier temperature control, 10 to 100 °C, for 10 mm cells | 634-0576 |
| Thermostatable single cell holder, thermostat required | 634-0059 |
| Micro cell holder for HITACHI micro cells | 634-0701 |
| Tandem cell holder for max. 3 cells on sample and reference position | 634-0062 |
| 5 cell changer, manual | 634-0058 |
| 4 position rectangular long path cell holder for 10, 20, 30, 40, 50, 100 mm | 634-0054 |
| Single rectangular long cell holder for 10, 20, 30, 40, 50, 100 mm | 634-0055 |
| Single long cell holder for cylindrical cells | 634-0056 |
| Glass filter holder | 634-0057 |
| Film holder | 634-0060 |
| Polarisation filter holder including polarisation filter | 634-0063 |
| Automatic 6 cell holder with Peltier temperature control | 634-0722 |
| Specular reflectance sample holder, 5°, relative | 634-0575 |
| Integrating sphere, 60 mm internal diameter | 634-0574 |
| Integrating sphere, 150 mm internal diameter | 634-0573 |

The U-3900 series are true dual beam spectrophotometers that have been specially developed for research and development.

These devices are highly sensitive, and have selectable bandwidth and a high level of performance. The U-3900H is particularly suited for measuring samples with high absorption, as the very low scattered light of the pre-monochromator ensures a high linearity across a broad absorption range. The sensitivity of the UV range can be increased by setting the scan rate in the UV range separately.

- Stable baseline in the wavelength range of 190 - 850 nm
- High accuracy and reproducibility thanks to the chopped true dual beam optics and photomultiplier detector
- UV Solutions® for U-3900 software provides flexibility for applications and reports
- Scan speed for the UV range can be set individually for increased sensitivity
- The USB interface for PC connection also allows modern laptops or notebooks to be used (that no longer have a RS232 interface)
- A notebook PC can be placed on the photometer to save space

Technical specifications

| Model | U-3900 | U-3900 H |
|-----------------------------|---|----------------------------|
| Light source | Tungsten and deuterium lamps | |
| Optical system | True dual beam with chopper and photomultiplier | |
| Monochromator | Single | Double |
| Sample space dimensions | Layer thickness: 100 mm W x D x H: 120 x 300 x 140 mm | |
| Wavelength range (nm) | 190 - 900 | |
| Scan speed (nm/min) | 1,5; 3; 15; 30; 60; 120; 300; 600; 1200; 1800; 2400 | |
| Spectral bandwidth (nm) | 0,1; 0,5; 1; 2; 4,5 | |
| Wavelength accuracy (nm) | ±0,1 (at 656,1 after wavelength calibration) | |
| Photometric measuring range | -3,8 to +3,8 A; 0 - 300% T | -5,5 to +5,5 A; 0 - 300% T |
| Photometric accuracy | ±0,002 A (0 - 0,5 A) ±0,003 A (0,5 - 1,0 A) ±0,006 A (1,0 - 2,0 A) ±0,3% T | |
| Stray light | 0,015% | 0,00025% |
| Interfaces | USB | |
| Power (V) | 100 - 240 | |
| W x D x H (mm) | 680 x 692 x 257 | |
| Weight (kg) | 45 (without PC) | |

Ordering information

| Description | Cat. No. |
|---|----------|
| U-3900 spectrophotometer, only EU plug version available | 634-0756 |
| U-3900H spectrophotometer, only EU plug version available | 634-0754 |

Fluorometers

Fluorescence spectrophotometer F-2700



F-2700 from HITACHI

Unique stand-alone scanning fluorometer offering quantum yield measurement options

The F-2700 fluorescence spectrophotometer surpasses the competition with its superior performance, functionality and flexibility. Within one instrument you can choose between stand-alone or PC-controlled function (together with optional software). Variable slit width and a wide dynamic range offers exceptional performance. Its ultra-high sensitivity (S/N: 800 RMS) combined with the scan speed of 3000 nm/min (stand-alone or 12 000 nm/min PC controlled) is unique in this class of instrument.

Above all those options the instrument offers a huge number of accessories including the option to measure the quantum yield of substances.

All this together in a very compact and space-saving unit.

Ordering information

| Description | Cat. No. |
|--|----------|
| Fluorescence spectrophotometer F-2700 with EU plug | 634-0067 |

Delivery information

Supplied with a xenon lamp, fluorescence cuvette, tool kit and user manual

Accessories

| Description | Cat. No. |
|---|----------|
| Thermostatable single cell holder | 634-0713 |
| Solid sample cell holder | 634-0075 |
| Polarisation filter accessories for the UV/VIS range from 260 to 700 nm | 634-0746 |
| Polarisation filter accessories for the VIS range from 380 to 800 nm | 634-0073 |
| Cut-off filter set | 634-0747 |

More accessories available on request.

- Excellent sensitivity (signal/noise ratio: 800 or better)
- Spectral bandwidth selectable from 2,5 to 20 nm
- Large dynamic linear measuring range
- Spectral correction (option under PC control) for comparable and true spectra
- Quantum yield accessory (option under PC control)
- Easy-to-use keypad
- A horizontal beam path allows small sample volumes (approx. 0,6 ml in standard cells) to be measured

Technical specifications

| | |
|--|--|
| Light source | 150 W xenon lamp |
| Sensitivity | Signal/noise ratio 800 or better (Raman of water: EX 350 nm, slit: 5 nm) |
| Display range | 0,000 - 9999 |
| Wavelength range (nm) | 220 - 730 (optional R928F 220 to 800) |
| Wavelength resolution (nm) | 2,5 |
| Spectral bandwidth (nm)(excitation and emission) | 2,5; 5; 10; 20 |
| Scan speed (nm/min) | 60, 300, 1500, 3000 (stand-alone) 60, 300, 1500, 3000, 12 000 (PC) |
| Interface | USB to PC |
| Display | Monochromatic 5,7" LCD |
| Sample holder | 10 mm rectangular cell holder |
| Power | 100 - 240 V AC; 50/60 Hz |
| W x D x H (mm) | 503 x 600 x 343 |
| Weight (kg) | 41 |

Ordering information

| Description | Cat. No. |
|--|----------|
| Fluorescence spectrophotometer, F-7000 including UV FL Solutions® software | 634-0743 |

Delivery information

(PC not supplied). System specifications are available on request.

Accessories

| Description | Cat. No. |
|---|----------|
| Sipper | 634-0736 |
| Single cell holder for sensitive measurements | 634-0737 |
| Thermostatable single cell holder | 634-0713 |
| Thermostatable single cell holder with stirring function | 634-0714 |
| 4x changer, manual | 634-0753 |
| 8x changer, manual | 634-0752 |
| Accessories for low cryogenic temperature measurements | 634-0735 |
| Sample holder for solids | 634-0075 |
| Cell holder for absorption measurements | 634-0748 |
| Accessories for measurements with micro volumes | 634-0734 |
| Polarisation filter accessories for the UV/VIS range from 260 to 700 nm | 634-0746 |
| Polarisation filter accessories for the VIS range from 380 to 800 nm | 634-0073 |
| Accessories for polarisation measurements UV/VIS, automatic | 634-0738 |
| Accessories for polarisation measurements VIS, automatic | 634-0739 |
| Cut-off filter set | 634-0747 |
| Microplate holder | 634-0740 |
| Intracellular cation measuring programme | 634-0742 |
| Report generator for FL Solutions® software | 634-0741 |

F-7000 from HITACHI

A superior device bearing the HITACHI name with excellent performance.

This fluorescence photometer offers comprehensive possibilities with outstanding scan speeds of 60 000 nm/min. It is especially suitable for 3-D measurements. The horizontal measuring beam also allows measurements to be taken for small volumes (from 0,6 ml in macro cuvettes).

The user-friendly FL Solutions® software is a powerful tool for analyses and provides a wide range of functions, including 3-D measurements and control of the optional associated microplate holder. Other additional applications are phosphorescence measurements, FRET, BRET and Ca²⁺ identification in cells.

- High sensitivity optical system for precise measurements of low sample volumes: Signal/noise 800 (RMS), 250 (peak to peak)
- Extremely rapid scan speed of 60 000 nm/min allows fast reactions to be recorded with 3-D measurements
- Large range of applications in the fields of material industry, pharmaceuticals and life science research
- Linear signal over up to 6 dimensions of the concentration



Technical specifications

| | |
|---|---|
| Light source | 150 W xenon lamp |
| Sensitivity | Signal/noise ratio 800 or better (RMS), signal/noise ratio 250 or better (peak to peak) |
| Monochromator | Stigmatic concave grating: 900 lines/mm |
| Detector | Photomultiplier |
| Display range of the display | -9,999 to +9,999 |
| Wavelength range (nm)(stimulation and emission) | 220 - 750 |
| Wavelength resolution (nm) | 1,0 (at 546,1) |
| Wavelength accuracy (nm) | ±1,0 |
| Spectral bandwidth (stimulation and emission) | 1; 2,5; 5; 10; 20 |
| Scan speed (nm/min) | 30, 60, 240, 1200, 2400, 12 000, 30 000, 60 000 |
| Interface | USB |
| PC software | FL Solutions® for F-7000 |
| Sample holder | 10 mm rectangular cell holder |
| Power (V) | 100 - 240 |
| W x D x H (mm) | 620 x 520 x 300 |
| Weight (kg) | 41 |

Turbidimeter equipment for measuring turbidity

Turbidimeter Turbiquant® 1100 IR/T - 1500 IR/T - 300 IR/T

The degree of turbidity is used as a measurement of water quality.

Cloudy liquid contains undissolved suspended matter (particles, bacteria etc.). These particles scatter incident light in many directions. The light scattered at a 90 degree angle is measured (nephelometry). The clearer the liquid, the lower the turbidity. The units of turbidity are measured as NTU = Nephelometric Turbidity Units. Drinking water has a maximum of 5 NTU, sewage more than 2000. The Turbiquant® devices perform every kind of turbidity measurement convincingly and deliver the results quickly, simply and accurately.

Technical specifications

| Model | Turbiquant® 1100 IR/1100 T | Turbiquant® 1500 IR/1500 T | Turbiquant® 3000 IR/3000 T |
|-------------------|---|--|--|
| Measurement units | NTU, FNU | NTU | NTU, FNU, FAU, EBC |
| Range | 0,01 - 1100 NTU | 0 - 1000 NTU | 0 - 10 000 NTU/0 - 10 000 FNU/0 - 10 000 FAU; 0 - 2450 EBC |
| Resolution | 0,01/0,1/1 depending on measurement range | | Selectable 0,1 - 0,0001 NTU 0,0001/0,001/0,01/0,1 depending on the measuring range |
| Accuracy | ±2% of measured value or ±0,1 NTU in the range 0 - 500 NTU ±3% of the measured value in the range 500 - 1100 NTU | ±2% of the measured value or ±0,1 NTU in the range 0,00-1000 NTU | ±2% of the measured value or ±0,1 NTU, the larger value of the two in the range 0,001000 NTU ±5% of the measured value in the range, 1000 - 4000 NTU ±10% of the measured value in the range 4000 - 10 000 NTU |
| Repeatability | - | <±1% of measured value or ±0,01 NTU, whichever is the larger | |
| Calibration | Automatic 1 - 3 point | | Automatic 1 - 4 point (to 1750 NTU) 1750 - 10 000 NTU selectable |
| Cuvettes | 25 x 45 mm | 28 x 70 mm or optional flow through cell | 28 x 70 mm or optional flow through or low pressure flow through cells |
| Volume (ml) | 15 | 25 | 25 |
| Interface | - | RS232 bi-directional | |
| Time/calendar | - | Integrated | |
| GLP functions | - | Control of calibration intervals, self-test | Control of calibration intervals, self-test, calibration and instrument configuration are access code protected |
| Power | 4x alkali manganese batteries (AAA) | Universal charger/plug | |
| Certification | CE | CE, UL, TÜV/GS | |
| Light source | For a model T tungsten halogen lamp as per USEPA method, for a model IR with infrared light source as per EN ISO 7027 | | |



Turbidimeter equipment for measuring turbidity

Turbidimeter Turbiquant® 1100 IR/T - 1500 IR/T - 300 IR/T

Turbiquant® 1100 IR/1100 T

- A portable, battery operated turbidity meter for on-site analysis
- Battery power for more than 5000 measurements
- Waterproof housing (IP67 compliance) with easy-to-use 5 button operation
- Rugged carrying case

Items supplied: Turbidimeter, 2 empty cells, quick guide, user manual, carrying case and 4 batteries.

Turbiquant® 1500 IR/1500 T

Ideal for routine analytical requirements in the laboratory

- Simple handling
- Automatically selects relevant resolution
- Optional pour through assembly for high sample throughput

Items supplied:

Turbidimeter, universal charger/plug, 3 empty cells, quick guide and user manual.

Turbiquant® 3000 IR and 3000 T for optimum flexibility

Same as the 1500 models plus the following features:

- Wider measuring range - up to 10 000 NTU
- 4 detectors give nephelometric and ratio measurements
- Transmission measurements over 40 FNU in line with EN ISO 7027
- Optional flow through assembly

Items supplied:

Turbidimeter, universal charger/plug, 3 empty cells, quick guide and user manual.

Ordering information

| Description | Cat. No. |
|--|--------------|
| Turbiquant® 1100 IR turbidimeter, portable plug type | 1.18324.0001 |
| Turbiquant® 1100 T turbidimeter, portable | 1.18325.0001 |
| Turbiquant® 1500 IR turbidimeter | 1.18330.0001 |
| Turbiquant® 1500 T turbidimeter | 1.18331.0001 |
| Turbiquant® 3000 IR turbidimeter | 1.18332.0001 |
| Turbiquant® 3000 T turbidimeter | 1.18333.0001 |

Accessories

| Description | Pk | Cat. No. |
|---|-------|--------------|
| Turbiquant® 1100 IR empty cuvettes | 3 | 1.18320.0001 |
| Turbiquant® 1500/3000 cuvettes | 3 | 1.18336.0001 |
| Lamp, LED for Turbiquant® 1500 IR | 1 | 1.18344.0001 |
| Turbiquant® 1500/3000 flow through cells | 1 | 1.18340.0001 |
| Turbiquant® 1500/3000 tungsten lamp | 1 | 1.18338.0001 |
| Lamp, LED for Turbiquant® 3000 IR | 1 | 1.18382.0001 |
| Turbiquant® 1500/3000 cuvette stands | 1 | 1.18339.0001 |
| Turbiquant® 3000 flow through cells (low pressure) | 1 | 1.18341.0001 |
| Printer cable for Turbiquant® 1500/3000 (serial) | 1 | 1.09759.0001 |
| PC cable for Turbiquant® 1500/3000 (for serial interface) | 1 kit | 1.14667.0001 |

Calibration standards

| Description | Pk | Cat. No. |
|--|----|--------------|
| Turbiquant® 1000 IR calibration standard set (0,02 - 10,0 - 100 - 1000 NTU) | 1 | 1.18327.0001 |
| Turbiquant® 1500 IR/1500 T calibration standard set (0,02 - 10,0 - 1000 NTU) | 1 | 1.18328.0001 |
| Turbiquant® 3000 IR calibration standard set (0,02 - 10,00 - 100,0 - 1750 NTU) | 1 | 1.18329.0001 |
| Turbiquant® 1100 IR/1100 T (0,02 - 10,0 - 1000 NTU) standard set | 3 | 1.18335.0001 |
| Turbiquant® 3000 IR calibration standard (10 000 NTU) | 1 | 1.18342.0001 |
| Turbiquant® 3000 T calibration standard (10 000 NTU) | 1 | 1.18343.0001 |
| Turbiquant® 1500 IR/1500 T calibration standard set (0,02 - 10,0 - 100,0 - 1750 NTU) | 1 | 1.18349.0001 |



Cuvettes

Absorption measurement macro standard cuvettes

Glass and quartz cuvettes from Hellma® for spectrophotometers

These cuvettes are manufactured from various types of glass. The most important criteria for the choice of the type of glass is the spectral range for which the cell is intended. Each cuvette is identified by an etched on colour symbol that indicates the spectral range.

The layer thickness plays a particularly important part in some photometric applications.

| Material | Layer thickness (mm) | Tolerance (± mm) |
|-----------------------|----------------------|------------------|
| Quartz glass | 0,01 to 0,05 | 0,003 |
| | 0,1 to 0,2 | 0,005 |
| | 0,5 to 20 | 0,01 |
| | 40 to 100 | 0,02 |
| Special optical glass | 0,1 to 20 | 0,01 |
| | 40 to 100 | 0,02 |
| Optical glass | 10 to 30 | 0,1 |
| | 40 to 100 | 0,2 |



Macro cuvette

Technical specifications

| Logo | Glassware type | Spectral range (nm) | Applications/spectrum |
|----------|--|---------------------|-----------------------|
| QS | High quality SUPRASIL® quartz glass | 200 to 2500 | UV |
| QH or UV | Natural crystalline quartz | 230 to 2500 | UV |
| QX | Synthetic quartz glass without OH absorption | 200 to 3500 | Near IR |
| OS | Special optical glass | 320 to 2500 | Visible |
| OG | Optical glass - not normally used in spectrophotometry | 360 to 2500 | Visible |
| BF | Chemical-resistant borosilicate glass | 330 to 2500 | Visible |

Ordering information

| Style | Window material | Layer thickness (mm) | Height (mm) | Volume (µl) | Cat. No. |
|---------|----------------------------|----------------------|-------------|-------------|----------|
| 6030-OG | Optical glass | <10 | 45 | 3500 | 634-9001 |
| | | 20 | 45 | 7000 | 634-9002 |
| | | <40 | 45 | 14 000 | 634-9003 |
| | | 50 | 45 | 17 500 | 634-9004 |
| 100-OS | Special optical glass | 1 | 45 | 350 | 634-9011 |
| | | 2 | 45 | 700 | 634-9012 |
| | | 5 | 45 | 1750 | 634-9013 |
| | | 10 | 45 | 3500 | 634-9014 |
| | | 20 | 45 | 7000 | 634-9015 |
| | | 40 | 45 | 14 000 | 634-9016 |
| | | 50 | 45 | 17 500 | 634-9017 |
| | | 100 | 45 | 35 000 | 634-9018 |
| 100-QS | SUPRASIL® quartz glass | 1 | 45 | 350 | 634-9021 |
| | | 2 | 45 | 700 | 634-9022 |
| | | 5 | 45 | 1750 | 634-9023 |
| | | 10 | 45 | 3500 | 634-9024 |
| | | 20 | 45 | 7000 | 634-9025 |
| | | 40 | 45 | 14 000 | 634-9026 |
| | | 50 | 45 | 17 500 | 634-9027 |
| | | 100 | 45 | 35 000 | 634-9028 |
| 100-QX | SUPRASIL® 300 quartz glass | 10 | 45 | 3500 | 634-9034 |
| 110-OS | Special optical glass | 1 | 52 | 350 | 634-9041 |
| | | 2 | 52 | 700 | 634-9042 |
| | | 5 | 46 | 1750 | 634-9043 |
| | | 10 | 46 | 3500 | 634-9044 |
| | | 20 | 46 | 7000 | 634-9045 |
| | | 40 | 46 | 14 000 | 634-9046 |
| | | 50 | 46 | 17 500 | 634-9047 |
| 110-QS | SUPRASIL® quartz glass | 1 | 52 | 350 | 634-9051 |
| | | 2 | 52 | 700 | 634-9052 |
| | | 5 | 46 | 1750 | 634-9053 |
| | | 10 | 46 | 3500 | 634-9054 |
| | | 20 | 46 | 7000 | 634-9055 |
| | | 40 | 46 | 14 000 | 634-9056 |
| | | 50 | 46 | 17 500 | 634-9057 |
| | | 100 | 46 | 35 000 | 634-9058 |
| 110-QX | SUPRASIL® 300 quartz glass | 10 | 46 | 3500 | 634-9064 |

semi-micro cuvette



Micro cuvette

Ultra-micro cuvette



Semi-micro standard cuvettes for absorption measurements

| Style | Window material | Layer thickness (mm) | Height (mm) | Volume (µl) | Cat. No. |
|------------|----------------------------|----------------------|-------------|-------------|----------|
| 6040-OG | Optical glass | 10 | 45 | 1400 | 634-9101 |
| 6040-UV | HERASIL® quartz glass | 10 | 45 | 1400 | 634-9223 |
| 104-OS | Special optical glass | 10 | 45 | 1400 | 634-9112 |
| | | 50 | 45 | 7000 | 634-9115 |
| 104-QS | SUPRASIL® quartz glass | 10 | 45 | 1400 | 634-9122 |
| | | 50 | 45 | 7000 | 634-9125 |
| 104B-OS* | Special optical glass | 10 | 45 | 1400 | 634-9140 |
| Semi-micro | SUPRASIL® 300 quartz glass | 10 | 45 | 1400 | 634-9132 |
| 104B-QS* | SUPRASIL® quartz glass | 10 | 45 | 1400 | 634-9150 |
| 114-OS** | Special optical glass | 10 | 46 | 1400 | 634-9148 |
| 114-QS** | SUPRASIL® quartz glass | 5 | 46 | 700 | 634-9156 |
| | | 10 | 46 | 1400 | 634-9158 |

* Black side walls and base

** Round PTFE stoppers

Micro cuvettes for absorption measurements

Micro cuvettes with PTFE lid, internal dimensions: width 2 mm from Hellma®

| Style | Window material | Layer thickness (mm) | Height (mm) | Volume (µl) | Cat. No. |
|--------------|------------------------|----------------------|-------------|-------------|----------|
| 104.002-OS | Special optical glass | 10 | 45 | 700 | 634-9079 |
| 104.002B-OS* | | 10 | 45 | 700 | 634-9214 |
| 104.002-QS | SUPRASIL® quartz glass | 10 | 45 | 700 | 634-9080 |
| 104.002B-QS* | | 10 | 45 | 700 | 634-9081 |

Ultra-micro cuvettes for absorption measurements

With PE stoppers.

| Style | Window material | Layer thickness (mm) | Height (mm) | Volume (µl) | Cat. No. |
|------------|------------------------|----------------------|-------------|-------------|----------|
| 105.200-QS | SUPRASIL® quartz glass | 10 | 45/15* | 180 | 634-9087 |
| | | 10 | 45/8,5** | 180 | 634-9088 |
| 105.201-QS | | 10 | 45/15* | 120 | 634-9089 |
| | | 10 | 45/8,5** | 100 | 634-9090 |
| 105.202-QS | | 10 | 45/15* | 70 | 634-9091 |
| | | 10 | 45/8,5** | 70 | 634-9092 |
| 105.203-QS | | 10 | 45/15* | 70 | 634-9328 |
| | | 10 | 45/8,5** | 70 | 634-9093 |

* Black side walls and base

** Round PTFE stoppers

Cuvettes

Disposable cuvettes VWR Collection



Disposable cuvettes for spectrophotometers

Disposable cuvettes made from optical PS with optimal transmittance from 340 to 900 nm or made from PMMA with an application range of 300 to 900 nm.

- Very low variation of absorbance values
- Layer thickness 10 mm, with 2 windows
- Packed in Styrofoam racks and wrapped in plastic film

Ordering information

| Style | Materials | Volume (ml) | Pk | Cat. No. |
|----------------------------|-----------|-------------|-----|----------|
| Macro | PS | 4 | 100 | 634-2500 |
| Semi-micro, VWR Collection | PS | 1,6 | 100 | 634-2501 |
| | PMMA | 4 | 100 | 634-2502 |
| Macro | PMMA | 1,6 | 100 | 634-2503 |

Disposable UV cuvettes Brand



Lids for UV micro cuvettes, PE

Ordering information

| Colour | Pk | Cat. No. |
|--------|-----|----------|
| Blue | 100 | 612-5690 |
| Yellow | 100 | 612-5691 |
| Green | 100 | 612-5692 |
| Orange | 100 | 612-5693 |

Disposable UV cuvettes

These disposable plastic UV cuvettes from BRAND can be used in many fields of application instead of expensive and sensitive cuvettes made from glass or quartz glass. Their high chemical resistance means they can be used with most polar solvents as well as with acids and alkalis. The UV micro cuvettes can be used for measurements from 220 nm and with small sample quantities (70 µl are enough for measurement). The cuvette layer thickness is 10 mm. Two different centre heights ensure they can be used in most commercially available photometers

- Specially developed for the photometric measuring of proteins, DNA and RNA in the UV range
- Ideally suited for measurements at 260 nm, 280 nm and in the visible wavelength range
- Colour circular closures for the micro cuvettes ensure that samples can be easily recognised and the cuvettes are reliably leak-proof at up to -20 °C

Ordering information

| Style | Volume | Pk | Cat. No. |
|--|--------------|-----|----------|
| UV cuvettes, semi-micro, 4,5 x 23 mm aperture | 1,5 - 3,0 ml | 100 | 612-5684 |
| UV cuvettes, macro, 10 x 35 mm aperture | 2,5 - 4,5 ml | 100 | 612-5685 |
| UV cuvettes, micro, z=8,5 mm, 2 x 3,5 mm aperture | 70 - 850 µl | 100 | 612-5686 |
| | 70 - 850 µl | 500 | 612-5687 |
| UV cuvettes, micro, z=8,5 mm, 2 x 3,5 mm aperture, individually packed | 70 - 850 µl | 100 | 612-2505 |
| UV cuvettes, micro, z=15 mm, 2 x 3,5 mm aperture | 70 - 550 µl | 100 | 612-5688 |
| | 70 - 550 µl | 500 | 612-5689 |
| UV cuvettes, micro, z=15 mm, 2 x 3,5 mm aperture, individually packed | 70 - 550 µl | 100 | 612-2506 |

UVette® 80


- Individually blister packed for contamination-free work
- DNA, RNAase and protein-free

UVette® Starter set

- 80 individually packed UVettes® and a universal adapter for 15 mm light beam height (including GeneQuant); can be converted to 8,5 mm

UVette® Routine pack

- 200 UVettes® bulk packed in a resealable box (as an economical solution)



UVette® disposable cuvettes for spectrophotometers

The plastic material of the disposable UVette® provides a transparency spectrum of 220 to 1600 nm. The cuvette has a funnel shaped tube bottom that prevents capillary effects and ensures that the liquid always remains in the centre of the measuring range; this is particularly important for micro volumes such as 50 µl. The design provides a 2 or 10 mm optical path length by simply rotating the UVette® by 90 degrees. The 2 mm layer thickness means that no dilution of highly concentrated samples is required in most cases.

- Compatible with most spectrophotometer models (with appropriate adapter)
- 4 optical windows and volume markings for 500 and 1000 µl
- No adapter required when used with the Eppendorf® BioPhotometer



Ordering information


| Description | Pk | Cat. No. |
|--|-------|----------|
| UVette® 80 | 80 | 634-1921 |
| UVette® Starter set | 1 kit | 634-1922 |
| UVette® Routine pack | 200 | 634-1923 |
| UVette® adapter for 8,5 mm light beam height (centre height) | 1 | 634-0022 |
| UVette® adapter for 10 mm light beam height (centre height) | 1 | 634-0023 |
| UVette® adapter for 15 mm light beam height (centre height) | 1 | 634-0024 |
| UVette® adapter for GeneQuant | 1 | 634-0025 |

Accessories - Cell stands, cleaning agents, software

Cell stands

PE, for 12 cuvettes with 10 mm layer thickness - Kartell

| Description | Cat. No. |
|-------------|----------|
| Cell stand | 634-8580 |



Cuvette cleaning agents

An alkaline liquid concentrate from Hellma®, which is typically diluted to 2% with water for cleaning quartz glass and glass cuvettes. It can also be used to clean other sensitive optical components made of glass, quartz, sapphire and porcelain. It is a complete replacement for chromo sulphuric acid. It removes even the most stubborn contaminants such as fats, waxes, and proteins whilst being considerably less hazardous and corrosive than many of the traditional alternatives.

| Description | Pk | Cat. No. |
|----------------|-----|----------|
| Hellmanex® III | 1 l | 634-0666 |



Quality assurance

Secondary standards for calibrating spectrophotometers

UV/ VIS spectrophotometry is one of the most frequently applied methods in chemical analysis. It is used in clinical chemistry, the pharmaceutical industry, scientific research and very often in quality assurance.

In recent years quality requirements, outlined by EN ISO 9000, Good Laboratory Practice (GLP), Good Manufacturing Practice (GMP) or the recommendations of pharmacopeias (EP, DAB, USP) have become more influential. Formal performance verification of UV/VIS spectrophotometers is now essential. The performance verification tests required by major Pharmacopeias for UV/VIS spectrophotometers check the spectral resolution, the wavelength accuracy, the photometric accuracy and for stray light.

A set of liquid Hellma® calibration standards enables you to check spectrophotometers with respect to wavelength and photometric accuracy, stray light and spectral resolution at wavelengths from 198 to 650 nm.

They also enable the checking of photometric accuracy in the visible range of the spectrum and the wavelength accuracy in the ultraviolet and visible range.



Solid filter set

The calibration set consists of three neutral density glass filters for checking absorbance accuracy, and a holmium oxide glass filter for checking wavelength accuracy.

The glass filters are mounted in precision aluminium frames and are designed for use with the standard 10 mm cell holder provided with spectrophotometers. The set, along with an empty mount, is supplied in a sturdy wooden case.

For identification purposes the set number is engraved on each filter mount. The absorbance values and/or peak position wavelength of every filter is quoted in the accompanying calibration certificate. A copy of the values, for laboratory use, can be found on the inside of the case lid.

Ordering information

| Description | Cat. No. |
|--|-------------------|
| Complete set for testing photometric and wavelength accuracy | 634-9691 |
| Holmium oxide glass filter to test wavelength accuracy | 634-9693 |
| Glass filter for testing photometric accuracy (nominal value of the absorption 0,25 A) | 634-9694 |
| Glass filter for testing photometric accuracy (nominal value of the absorption 0,5 A) | 634-9695 |
| Glass filter for testing photometric accuracy (nominal value of the absorption 1 A) | 634-9696 |
| Didymium glass filter for testing photometric accuracy and wavelength accuracy | HELL666-000F7-323 |

Liquid filter set

Consisting of eight liquid filters which enable the complete checking of the spectrophotometer according to European Pharmacopeia. They are designed for use with a standard 10 mm cell holder provided with spectrophotometers. The set is supplied in a sturdy wooden case. For identification purposes an ID number engraved on each filter. The absorbance values and/or peak position wavelengths of every filter are quoted in the accompanying calibration certificate.

Ordering information

| Description | Pk | Cat. No. |
|--|-------|-------------------|
| Set for checking spectrophotometer according to Ph. Eur. | 1 kit | HELL667-003-UV-40 |
| Set for testing stray light according to Ph. Eur. | 1 | HELL667-100-UV-40 |
| Set for testing the resolution according to Ph. Eur. | 1 | HELL667-200-UV-40 |
| Set for testing the photometric accuracy according to Ph. Eur. | 1 | HELL667-305-UV-40 |

Absolutely seamless

Consistent quality assurances convert your measurements into safe and verifiable analysis results. The Spectroquant® NOVA and Pharo photometers ensure your work is GLP compliant. The building blocks of the AQS concept are tried and tested media for your Internal Quality Control (IQC), as recommended in the data sheet ATV-DWVK A 704 of the German Association for Water, Wastewater and Waste (ATV-DWVK).



The choice is yours:

AQS-1 mode for monitoring the photometer

The photometer is monitored in AQS-1 mode using certified colour standards (Spectroquant® PhotoCheck) and also for the Spectroquant® Pharo spectrophotometers using CertiPUR® UV/VIS standards and Spectroquant® Photocheck.

Spectroquant® PhotoCheck is a complete set of long-term stable dye solutions. The traceability of this test equipment on international standards is guaranteed by testing in a reference spectrophotometer. This is monitored using primary standards (NIST Standards). This means that Spectroquant® PhotoCheck is traceable and therefore suitable for monitoring test equipment in accordance with DIN ISO 9001 or 14001. All results can ultimately be transferred to a printer or PC for documentation purposes.

Ordering information

| Description | Pk | Cat. No. |
|--------------------------|-------|--------------|
| Spectroquant® PhotoCheck | 1 kit | 1.14693.0001 |

CertiPUR® UV/VIS standards

CertiPUR® standards can be used to check your spectrophotometer is functioning uniformly and correctly. The following parameters can be tested with the CertiPUR® solutions in accordance with Ph. Eur.:

- Absorbancy
- Scattered light behaviour
- Wavelength accuracy

These regular checks are required in accordance with GLP, GMP, USP and DIN 9001 or EN 45001.

Ordering information

| Description | Volume | Cat. No. |
|--------------------|---|--------------|
| UV/VIS Standard 1 | Potassium dichromate solution for the absorption as per DAB and Ph. Eur. 2 x 10 ml $K_2Cr_2O_7$ 60,06 mg/l in H_2SO_4 0,01 N and 6 x 10 ml H_2SO_4 0,01 N | 1.08160.0001 |
| UV/VIS Standard 1A | Potassium dichromate solution for the absorption at 430 nm as per DAB and Ph. Eur. 2 x 10 ml $K_2Cr_2O_7$ 600,06 mg/l in H_2SO_4 0,01 N and 6 x 10 ml H_2SO_4 0,01 N | 1.04660.0001 |
| UV/VIS Standard 2 | Sodium nitrite solution for scattered light testing as per DAB and Ph. Eur. 3 x 10 ml $NaNO_2$ 50 g/l in H_2O | 1.08161.0001 |
| UV/VIS Standard 3 | Sodium iodide solution for scattered light testing as per DAB and Ph. Eur. 3 x 10 ml NaI 10 g/l in H_2O | 1.08163.0001 |
| UV/VIS Standard 6 | Holmium oxide solution reference material for the wavelength as per DAB and Ph. Eur. 3 x 10 ml HO_2O_3 40 g/l in $HClO_4$ (10% v/v) | 1.08166.0001 |



Quality assurance

Spectroquant® quality assurance

AQS-2 mode for monitoring the complete system

AQS mode monitors the complete system comprehensively using certified multi-parameter standards.

Spectroquant® CombiCheck

Spectroquant® CombiCheck solutions are ready-to-use multi-parameter standards for use with Spectroquant® test kits. Each pack contains a standard solution and an additional solution. Both solutions can be used immediately for quality assurance without requiring dilution.

The standard solution is used to prove the correctness of the results of the complete Spectroquant® system (photometer, reagents, analytical procedure and operation).

The additional solution is used to test sample-related influences by measuring the recovery rate and identifies the sample-related preparation.

The maximum number of identifications depend on the test kit used when applying the standard solution.

This function has additional device driven support for the Spectroquant® Pharo spectrophotometers.

Technical specifications

| Spectroquant® CombiCheck 10 | |
|-----------------------------|--|
| Ammonium | 4,00 ±0,30 mg/l NH ₄ ⁺ |
| Chloride | 25 ±6 mg/l Cl ⁻ |
| COD | 80 ±12 mg/l COD |
| Nitrate | 2,5 ±0,25 mg/l NO ₃ ⁻ |
| Phosphate | 0,80 ±0,08 mg/l PO ₄ ⁻ |
| Sulphate | 100 ±15 mg/l SO ₄ ²⁻ |
| Spectroquant® CombiCheck 20 | |
| Ammonium | 12,0 ±1,0 mg/l NH ₄ ⁺ |
| Chloride | 60 ±10 mg/l Cl ⁻ |
| COD | 750 ±75 mg/l COD |
| Nitrate | 9,0 ±0,9 mg/l NO ₃ ⁻ |
| Phosphate | 8,0 ±0,7 mg/l PO ₄ ⁻ |
| Sulphate | 500 ±75 mg/l SO ₄ ²⁻ |
| Spectroquant® CombiCheck 30 | |
| Cadmium | 0,500 ±0,060 mg/l Cd |
| Iron | 1,00 ±0,15 mg/l Fe |
| Copper | 2,00 ±0,20 mg/l Cu |
| Manganese | 1,00 ±0,15 mg/l Mn |
| Spectroquant® CombiCheck 40 | |
| Aluminium | 0,75 ±0,08 mg/l Al |
| Lead | 2,00 ±0,20 mg/l Pb |
| Nickel | 2,00 ±0,20 mg/l Ni |
| Zinc | 2,00 ±0,40 mg/l Zn |

| Spectroquant® CombiCheck 50 | |
|-----------------------------|--|
| Ammonium | 1,00 ±0,10 mg/l NH ₄ ⁺ |
| COD | 20,0 ±4,0 mg/l COD |
| Nitrogen | 5,0 ±0,70 mg/l N |
| Spectroquant® CombiCheck 60 | |
| Chloride | 125 ±13 mg/l Cl ⁻ |
| COD | 250 ±25 mg/l COD |
| Spectroquant® CombiCheck 70 | |
| Ammonium | 50 ±5 mg/l NH ₄ ⁺ |
| COD | 5000 ±400 mg/l COD |
| Nitrogen | 50 ±7mg/l N |
| Spectroquant® CombiCheck 80 | |
| COD | 1500 ±150 mg/l CSB |
| Nitrate | 25,0 ±2,5 mg/l NO ₃ ⁻ |
| Phosphate | 15,0 ±1,0 mg/l PO ₄ ⁻ |

Ordering information

| Description | Cat. No. |
|-----------------------------|--------------|
| Spectroquant® CombiCheck 10 | 1.14676.0001 |
| Spectroquant® CombiCheck 20 | 1.14675.0001 |
| Spectroquant® CombiCheck 30 | 1.14677.0001 |
| Spectroquant® CombiCheck 40 | 1.14692.0001 |
| Spectroquant® CombiCheck 50 | 1.14695.0001 |
| Spectroquant® CombiCheck 60 | 1.14696.0001 |
| Spectroquant® CombiCheck 70 | 1.14689.0001 |
| Spectroquant® CombiCheck 80 | 1.14738.0001 |



(AQS-3 marking)

AQS-1 + AQS-2:

The photometer indicates when the next AQS check is due. The intervals can be freely selected depending on time and also, for AQS-2, depending on the number of measurements.

AQS-3 matrix check

In addition to testing the complete system, measurement errors that can be caused by possible interferences in the sample also have to be identified.


Standard addition (see also Spectroquant® CombiCheck R-2 additional solutions) or dilution can be used to identify measurement errors that are caused by matrix influences. Recovery can be used to analyse these interferences and they can be removed using suitable measures e.g. a sample preparation.

Standard solutions

You can, of course, make use of our comprehensive range of ready-to-use standard solutions. Standards can be produced in the required concentration for almost every parameter using simple dilution.

Ordering information for standard solution (*traceable on NIST) 1000 mg/l


| Description | Cat. No. |
|-------------|--------------|
| Ammonium | 1.19812.0500 |
| Chloride | 1.19897.0500 |
| COD | 1.11769.0100 |
| Nitrate | 1.19811.0500 |
| Phosphate | 1.19898.0500 |
| Sulphate | 1.19813.0500 |



Water and Food Analytics

Test Kits, Instruments and Accessories

Quote W285102 for your own copy of this new catalogue



Technical services

We would be delighted to help you with the calibration, maintenance and repair of your testing equipment. Please contact your local sales office for further information.



Austria

VWR International GmbH
Graumanngasse 7
1150 Wien
Tel.: 01 97 002 0
Fax: 01 97 002 600
E-mail: info@at.vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: 016 385 011
Fax: 016 385 385
E-mail:
customerservice@be.vwr.com

Denmark

VWR - Bie & Berntsen
Transformervej 8
2730 Herlev
Tel.: 43 86 87 88
Fax: 43 86 87 90
E-mail: info@dk.vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: 09 80 45 51
Fax: 09 80 45 52 00
E-mail: info@fi.vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30 (0,15 EUR TTC/min)
Fax: 0 825 02 30 35 (0,15 EUR TTC/min)
E-mail: info@fr.vwr.com

Germany

VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Tel.: 0180 570 20 00*
Fax: 0180 570 22 22*
E-mail: info@de.vwr.com
*0,14 €/Min. aus d. dt. Festnetz,
Mobilfunk max. 0,42 €/Min.

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: (52) 521-130
Fax: (52) 470-069
E-mail: info@hu.vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20163 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com
info@internationalpbi.it

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

Labart Sp. z o.o.
A VWR International Company
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: labart@pl.vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

China

VWR International China Co., Ltd
Suite 3B02, Qilai Building, No. 889
Yishan Road
Shanghai 200233, China
Tel.: +86- 21 521 388 22
Fax: +86- 21 521 33 933
E-mail: sales_china@vwr.com

India

VWR Lab Products Pvt. Ltd
2nd Floor, Front Wing, 135/12, Brigade
Towers
Brigade Road
Bangaluru 560025 India
Tel.: +91-2522-647911/922
(Mumbai)
Tel.: +91-80-41117125/26 (Bangalore)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

GO TO VWR.COM FOR
THE LATEST NEWS,
SPECIAL OFFERS AND
DETAILS OF YOUR LOCAL
VWR DISTRIBUTOR