



ACE® HPLC / UHPLC Columns Specification Table

HPLC and UHPLC Columns

PHASE	USP LISTING	FUNCTIONAL GROUP	ENDCAPPED	PARTICLE SIZE (µm)*	PORE SIZE (Å)	SURFACE AREA (m ² /g)	CARBON LOAD (%)	PH RANGE	100% AQ compatible
ACE® Traditional Chemistries									
C18	L1	Octadecyl	Yes	1.7, 2, 3, 5, 10	100	300	15.5	2 – 8	-
C18-HL	L1	Octadecyl	Yes	3, 5, 10, 15	90	400	20.0	2 – 8	-
C8	L7	Octyl	Yes	2, 3, 5, 10	100	300	9.0	2 – 8	-
C4	L26	Butyl	Yes	2, 3, 5, 10	100	300	5.5	2 – 8	-
Ph	L11	Phenyl	Yes	2, 3, 5, 10	100	300	9.5	2 – 8	-
CN	L10	Cyano	Yes	2, 3, 5, 10	100	300	5.5	2 – 7	-
AQ	L1	Proprietary	Yes	2, 3, 5, 10	100	300	14.0	2 – 8	Yes
SIL	L3	Unbonded	No	2, 3, 5, 10	100	300	N/A	2 – 7	-
ACE NH ₂	L8	Proprietary aminopropyl	Proprietary	1.7, 3, 5	100	300	3.5	2 – 7	Yes
ACE® Novel Chemistries									
SuperC18	L1	Octadecyl encapsulated	Encapsulated	1.7, 2, 3, 5, 10	90	400	14.8	1.5 – 11.5	-
C18-AR	L1	Octadecyl with integral phenyl	Yes	1.7, 2, 3, 5, 10	100	300	15.5	2 – 8	Yes
C18-PFP	L1	Octadecyl with integral PFP	Yes	1.7, 2, 3, 5, 10	100	300	14.3	2 – 8	Yes
C18-Amide	L1 / L60	Polar embedded amide	Yes	1.7, 2, 3, 5, 10	100	300	16.4	2 – 8	Yes
CN-ES	L10	CN with extended alkyl spacer	Yes	1.7, 2, 3, 5, 10	100	300	12.6	2 – 8	Yes
Solid Core Technology ACE® Phases									
UltraCore SuperC18	L1	Octadecyl encapsulated	Encapsulated	2.5 5	95	130 100	7.0 5.4	1.5 – 11.0	-
UltraCore SuperPhenylHexyl	L11	Phenyl Hexyl encapsulated	Encapsulated	2.5 5	95	130 100	4.6 3.6	1.5 – 11.0	-
Large Molecule Wide Pore ACE® Phases									
C18-300	L1	Octadecyl	Yes	3, 5, 10	300	100	9.0	2 – 8	Yes
C8-300	L7	Octyl	Yes	3, 5, 10	300	100	5.0	2 – 8	Yes
C4-300	L26	Butyl	Yes	3, 5, 10	300	100	2.6	2 – 8	Yes
CN-300	L10	Cyano	Yes	3, 5, 10	300	100	2.6	2 – 7	Yes
Ph-300	L11	Phenyl	Yes	3, 5, 10	300	100	5.3	2 – 8	Yes
ACE® HILIC Phases									
ACE HILIC-A	L3	Proprietary SIL	No	1.7, 3, 5	100	300	N/A	2 – 7	-
ACE HILIC-B	L8	Proprietary aminopropyl	No	1.7, 3, 5	100	300	4.0	2 – 7	-
ACE HILIC-N	pending	Proprietary polyhydroxy	No	1.7, 3, 5	100	300	7.0	2 – 7	-
ACE® Method Development Kits (MDKs)									
ACE Advanced MDK	3 column kit	Contains C18, C18-AR, C18-PFP columns	Yes	1.7, 2, 3, 5	Data for these categories are as above for the individual column chemistries				
ACE Extended MDK	3 column kit	Contains SuperC18, C18-Amide, CN-ES columns	Yes	1.7, 2, 3, 5	Data for these categories are as above for the individual column chemistries				
ACE UltraCore MDK	2 column kit	Contains SuperC18, SuperPhenylHexyl columns	Yes	2.5, 5	Data for these categories are as above for the individual column chemistries				
ACE Biomolecule MDK	3 column kit	Contains C18-300, C4-300, Ph-300 columns	Yes	3, 5	Data for these categories are as above for the individual column chemistries				
ACE HILIC MDK	3 column kit	Contains HILIC-A, HILIC-B, HILIC-N columns	No	1.7, 3, 5	Data for these categories are as above for the individual column chemistries				